

INSTRUCTIONS FOR USE

CIC/MIC/ITC/ITE



## Thank you

Thank you for choosing our product as your means to better hearing. To support your efforts, we have put great care and attention into making sure that your new hearing instruments are of the highest quality and that they are easy to use and maintain.

We recommend that you read this manual carefully to achieve the maximum benefit of your new hearing instruments.

To meet your needs, various functions may be configured to your hearing instruments. The actual configuration of your instruments are marked (ticked off) by your Hearing Care Professional in relevant sections of this manual. Also see last page for a complete overview of all functions configured to your hearing instrument(s).

The hearing instrument amplification is uniquely adjusted and optimized to your personal hearing capabilities by the prescription performed by your Hearing Care Professional.

If you have any questions to the use or maintenance of hearing instruments, please contact your Hearing care professional.

## Indications for Use

The hearing instruments are intended to amplify and transmit sound to the ear and hereby compensate for mild-to-severe-to-profound impaired hearing.

## Contents


Instrument View CIC/MIC	7a
Instrument View ITC	7b
Instrument View ITC Power	7c
Instrument View ITE	7d
Warnings	9
Activating the hearing instrument for the first time	14
Changing batteries	15
Turning your instruments ON and OFF	17
Left/Right ear marking	18
Inserting your hearing instruments	19
Push-button (Optional)	20
Programs (Optional)	22
Volume Control by Push-button (Optional for ITC* only)	24
Mute (Optional)	25
Volume Control by Wheel (Optional for ITE and ITC Power)	26
Telecoil (Optional)	29

Auto Phone (Optional)	30
Caring for your hearing instruments	32
Microphone Protection System (Optional for CIC only)	38
Avoiding heat, humidity and chemicals	40
Using your hearing instrument	41
Common problems and their solutions	46
Wireless accessories	48
International Warranty	49
Warranty Certificate	50
Mobile Phone	51
Technical information	53
Settings overview for your hearing instrument	56


### IMPORTANT NOTICE

Please familiarize yourself with the entire contents of this booklet before using your hearing instruments. It contains instructions and important information about the use and handling of your hearing instruments and batteries.

## Instrument View CIC/MIC



Battery size: 10




*Instrument is not wireless*

### **IMPORTANT NOTICE**

Open the battery door fully to allow air to circulate whenever you are not using your hearing instrument, especially at night!

## **Instrument View ITC**




Battery size: 312


### **IMPORTANT NOTICE**

Open the battery door fully to allow air to circulate whenever you are not using your hearing instrument, especially at night!

## **Instrument View ITC Power**



Battery size: 312




*Instrument is not wireless*

### **IMPORTANT NOTICE**

Open the battery door fully to allow air to circulate whenever you are not using your hearing instrument, especially at night!

## **Instrument View ITE**



Battery size: 312

## Warnings

You should familiarize yourself fully with the following general warnings and the entire contents of this booklet before using your hearing instrument.

### Usage of hearing instruments

- Hearing instruments should be used only as directed and adjusted by your Hearing Care Professional. Misuse can result in sudden and permanent hearing loss.
- Never allow others to wear your hearing instrument as incorrect/wrongful usage could cause permanent damage to their hearing.

### Choking Hazards

- Hearing instruments, their parts, and batteries are not toys and should be kept out of reach of children and anyone who might swallow these items or otherwise cause injury to themselves.
- Never change the battery or adjust the controls of the hearing instruments in front of infants, small children or people with learning difficulties.
- Discard batteries carefully in a place where infants, small children or people with learning difficulties cannot reach them.

### IMPORTANT NOTICE

Open the battery door fully to allow air to circulate whenever you are not using your hearing instrument, especially at night!



- Batteries have occasionally been mistaken for pills. Therefore check your medicine carefully before swallowing any pills.
- Never put your hearing instrument or batteries in your mouth for any reason, as they are slippery and could be swallowed by accident.
- Most hearing instruments can be supplied with a tamper-resistant battery compartment upon request. This is strongly recommended for infants, small children, and people with learning difficulties.

*If a battery or hearing instrument is swallowed, see a doctor immediately.*

#### **Battery use**

- Always use batteries recommended by your Hearing Care Professional. Batteries of low quality may leak and cause bodily harm.
- Never attempt to recharge your batteries. They may explode and cause serious injury.
- Never dispose of batteries by burning them. There is a risk that they will explode and cause serious injury.

#### **Dysfunction in hearing instruments**

- Hearing instruments may stop functioning, for instance if the batteries have expired or if the tubing is blocked by moisture or earwax. You should be aware of this possibility, in particular when you are in traffic or otherwise dependent on warning sounds.

#### **Active Implants**

- Oticon shows caution and advise to follow the guideline as recommended by manufacturers of implantable defibrillators and pacemakers regarding use of mobile phones.
- If you wear an active implant, then keep the hearing instruments more than 15 cm away from the implant e.g. do not carry them in a breast pocket.  
If you have an active brain implant, please contact the manufacturer of your implantable device for information about the risk of disturbance.
- The included case for instrument storage has a built-in magnet. If you have an implantable device, such as pacemakers or defibrillators, the instrument case should not be carried in a chest pocket or near to the chest.

### **Explosives**

- The power source in your hearing instrument has insufficient energy to cause fire in normal conditions of use. The hearing instrument has not been tested for compliance with international standards concerning explosive atmospheres. It is recommended not to use your hearing instruments in areas where there is a danger of explosions.

### **X-ray, CT, MR, PET scanning and electrotherapy**

- Remove your hearing instrument for example during X-ray, CT / MR / PET scanning electrotherapy or surgery as your hearing instruments may be damaged when exposed to strong fields.

### **Possible side effects**

- Hearing instruments or ear moulds may cause an accelerated accumulation of ear wax.
- The otherwise non-allergenic materials used in hearing instruments may in rare cases cause a skin irritations or any other unusual condition.

*Please, seek consultation with a physician if these conditions occur.*

### **Interference**

- Your hearing instrument has been thoroughly tested for interference, according to the most stringent international standards. However, interference with your hearing instrument and other devices may occur, for example for some mobile telephones, citizens band devices and shop alarm systems. If this occur increase the distance between the hearing instrument and the device.

### **Wax protection system**

- If the instrument is equipped with a ProWax wax protection system, follow the directions specified in the section "How to change Wax filter".
- Always use the same type of wax filter like that which was originally supplied with the instrument.

*If you are in any doubt about the use or replacement of your wax protection system, contact your Hearing Care Professional.*

## Activating the hearing instrument for the first time

A hearing instrument is a miniature electronic device that runs on special batteries. To activate the hearing instrument, you must insert a fresh battery in the battery compartment.

The hearing instrument will start up when the battery door is closed. It will take a few seconds until sound is produced. A jingle is played while the start up process is running. Put on the hearing instruments upon closing the battery door.

If the hearing instrument is held in the hand while activated, it may produce a whistling sound. This will stop when properly fitted in the ear.

## Changing batteries


The correct battery size appears in the instrument overview in the first section of this document. Batteries perform differently. Ask your Hearing Care Professional for help choosing a high quality battery for your hearing instrument.

A worn out battery should be removed immediately and disposed according to local regulations.

You will hear 2 beeps when it is time to change the battery. This pre-warning will be repeated in moderate intervals until the battery runs out. The battery may run out shortly after a pre-warning or within the next couple of hours depending on battery capabilities and individual hearing losses.

If the battery runs out the instrument turns off. When that happens, 4 beeps are played to indicate that the instrument is no longer operating.


Replace the battery.



To replace the battery, follow these instructions:

- Open the battery door with your fingernail and remove the old battery. To assist with the handling of the batteries, the magnet at the end of the MultiTool can be used to help you.
- Remove the sticky label from the + side of the new battery.
- The new battery may be inserted by using the magnetic end of the MultiTool. Ensure that the positive (+) marking faces the positive (+) sign printed on the battery door. Close the battery door.

If any moisture is on the battery surface, it should be wiped off before use. When you change the battery, it will take a few seconds before the battery works.




## Turning your instruments ON and OFF

Turn **ON** your hearing instrument by closing the battery door completely with the battery in place.

When the battery door is closed the instrument will start up by playing a jingle. This indicates that the battery is working and the instrument is operating.

Turn **OFF** your hearing instrument by opening the battery door with your finger nail.

To preserve the battery, make sure your hearing instrument is switched off when you are not wearing it.




## Left/Right ear marking

Hearing instruments are fitted to the uniqueness of each ear, which means if you have two instruments then your left hearing instrument is shaped and programmed differently from your right. That is why it is important to distinguish between the left hearing instrument and the right.

In order to easily identify whether it is the left ear or right ear instrument, colour markings may be applied.

Blue text or a blue dot marks the **LEFT** instrument.  
Red text or a red dot marks the **RIGHT** instrument.




## Inserting your hearing instruments

Place the tip of the hearing instrument in your ear canal.

Gently pull your ear and push the hearing instrument into the ear canal, twisting slightly if necessary, following the natural contour of the ear canal.

*DO NOT use the battery door as a handle to insert or remove your instruments. It is not designed for this purpose.*



## Push-button (Optional)


Your hearing instrument may have a push-button. (If your hearing instrument does not have a push-button, this and the following 3 chapters are not relevant to you)

If you have one instrument, the push-button can be used for switching between different programs.

If you have two instruments, the push-button can be used for both adjusting the volume and switching between programs.\*

Finally, the push-button can be used to mute the instrument.

*\*Does not apply to CIC, MIC and ITC Power instruments*



The push-button on your hearing instrument can be programmed to one of 3 functionalities.

Left	Right	
		<b>Program shift:</b> Short press to change program.
		<b>Volume control change:</b> *Short press to adjust the volume.
		<b>Volume control and Program change*</b> Short press to adjust the volume and long push (about 2 seconds) to shift the program.

\*Only if you have two wireless instruments.

### Using 2 instruments

If you have two instruments the push-button on one instrument will change the status on both hearing instruments. This means, when you adjust the program or volume for one instrument, it automatically adjusts the other. (This feature does not apply to CIC, MIC and ITC Power instruments.)





*This is programmed by your Hearing Care Professional*

## Programs (Optional)

Your hearing instrument can have up to 4 different programs.

The number and type of programs available in your hearing instrument is programmed by your Hearing Care Professional.

When switching between the different programs, your hearing instrument will beep. The number of beeps will indicate which program you are using.

	One beep, when you switch to program 1
	Two beeps, when you switch to program 2
	Three beeps, when you switch to program 3
	Four beeps, when you switch to program 4

Press the push-button to switch forwards in the program cycle, e.g. from P1-P2.

If you have two instruments, then:

The **RIGHT** hearing instrument switches the program cycle forward, e.g. P1-P2.  
The **LEFT** hearing instrument switches program cycle, backwards e.g. P2-P1.

An overview and description of the hearing programs available in your hearing instrument is placed on the last page of this manual.

## Volume Control by Push-button (Optional for ITC\* only)

Only when using two instruments a Volume Control by Push-button can be activated. The volume control allows you to adjust the loudness of your hearing device given specific listening situations.

By pressing the push-button on one instrument, the volume is adjusted in both instruments simultaneously.

To **increase volume**, give a short press to the push-button on the RIGHT ear instrument.  
To **decrease volume**, give a short press to the push-button on the LEFT ear instrument.

You hear a beep when you increase or decrease the volume.

When the hearing instrument is turned on, it starts up at a preferred volume level. When adjusting the volume control to the preferred level, this is indicated by a beep.


*\*Does not apply to ITC Power instruments.*

## Mute (Optional)

Use the mute function if you need to silence the instrument while you have the instrument on.

Push the button for a minimum of 3 seconds to mute the instrument.

To reactivate the instrument, push the button briefly again.




### IMPORTANT NOTICE

Do not use the mute function as a switch off as the hearing instrument still draws current from the battery in this mode.



## Volume Control by Wheel (Optional for ITE and ITC Power)

Your hearing instrument may have a Volume Control wheel which always allows you to adjust the volume of each instrument. The Volume Control wheel can be turned with the tip of a finger and can be rotated infinitely both ways.



The Volume Control allows you to adjust the volume in specific listening situations to the level you feel comfortable with.

To **increase volume**, turn the wheel forward.

To **decrease volume**, turn the wheel towards the back of your head.

When the hearing instrument is turned on, it starts up at a preferred volume level. When adjusting the volume control to the preferred level, this is indicated by a beep.

## Using 2 instruments

If you have two instruments the Volume Control Wheel on one instrument can be programmed to adjust the volume on both hearing instruments (Binaural Coordination). This means, when you adjust the volume for one instrument, it automatically adjusts the other.\*

Binaural Coordination enabled.

*\*Not available for ITC Power instruments.*

## Telecoil (Optional)



Your hearing instrument may have a Telecoil.

The telecoil is a receiver for audio signals transmitted from an inductive loop installation. The telecoil is meant for telephone\* conversations and for loop systems in places such as theaters, churches or lecture rooms.



The telecoil is activated by the push-button. When activated, you will hear a certain number of beeps corresponding to the telecoil program.

See the “Hearing programs description” on last page for where the telecoil program is placed.

*\*Requires a telephone with built in tele-loop. Tele-loop is a common standard for assisted listening and transmits the audio to the telecoil in the hearing instrument.*


## Auto Phone (Optional)



The hearing instrument may have a built-in Auto Phone functionality. When the hearing instrument is close to a telephone receiver, the Auto Phone will activate a Phone Program. When the Phone Program is activated you will hear a number of beeps.

When you end your telephone conversation, the hearing instrument will automatically return to the previous program.

Not all telephones can activate the Auto Phone. The telephone receiver must have a special magnet fitted. Please follow the separate instructions for mounting the magnet, which is available from your Hearing Care Professional.



### WARNING

If you use an Autophone magnet:

- Keep magnets out of reach of children and pets. If a magnet is swallowed, see a doctor immediately.
- Do not wear the magnet in a breast pocket and always keep it 30 cm away from active implanatable devices such as pacemakers, defibrillators etc. Preferably, use the telephone with magnet on the opposite side than a pacemaker or other active implantants.
- Keep the magnet 30cm away from credit cards and other magnetically sensitive devices.


## Caring for your hearing instruments

Your ear canal produces ear wax (cerumen) which may clog up the sound outlet or the ventilation opening of your hearing instrument. Therefore your hearing instrument may contain a filter or another device to prevent wax accumulation. Consult your Hearing Care Professional for instructions about the specific wax protection system in your hearing instruments.

When handling a hearing instrument, keep it over a soft surface to avoid damage if you drop it.

## Cleaning the instrument

A range of dedicated cleaning tools is prepared for the special purpose of cleaning the instrument. These should be used to ensure the best care and performance of your hearing instruments.





The MultiTool contains a brush and a wire loop for the special purpose of cleaning wax from the hearing instrument.

### **IMPORTANT NOTICE**

The hearing instrument itself must never be washed or immersed in water or other liquids!

### Here is how . . .

- Clean away any traces of wax from around the ventilation opening and sound outlet with the brush. Hardened sediments can be removed by the wire loop in the MultiTool.
- If microphone openings are clogged by dust or dirt, this will affect the performance of the instrument. Carefully brush away visible signs of debris from the openings while twisting the brush.
- Clean the ventilation opening by pressing the brush through the hole while twisting it slightly.
- If the ventilation opening is very tiny a special tool may be required for cleansing. Your Hearing Care Professional can advise.



- Use the cloth to wipe clean the surface of the instrument.
- If the white wax protection system (called ProWax) placed in the sound outlet seems filled up and blocked by ear wax, it must be replaced.

Replace the MultiTool brush when necessary. Just pull it out of the handle and insert a new fresh brush. Press the new one firmly into the handle.





### Wax Protection System

- If the instrument is equipped with a ProWax wax protection system, follow the directions specified in the section "How to replace ProWax".

*If you are in any doubt about the use or replacement of your wax protection system, contact your Hearing Care Professional.*

## How to replace ProWax


1. Remove the tool from the shell. The tool has two arms, one with the new filter and one with the removal tool.
2. Push the removal tool into the existing waxfilter and remove the used filter from the hearing instrument.
3. Push the new filter into the hearing instrument.
4. Throw out the used tool.




## Microphone Protection System (Optional for CIC only)

### How to change T-Cap filter


1. Remove the tool from its packaging. The tool has two ends, one for removal of the T-Cap (A) and one with the new T-Cap (B).
2. Using the removal end, place it under the T-Cap edge, and push it under the top edge of the used T-Cap.
3. Lift the used T-Cap up and out
4. Using the other end of the tool, insert the new T-Cap into the opening and slightly twist the handle so that the new T-Cap is free.
5. Throw out the handle.




①




②



③



④



## Avoiding heat, humidity and chemicals

Your hearing instrument must never be exposed to extreme heat e.g. left inside a parked car in the sun. They must never be exposed to a lot of moisture e.g. steam baths, showers or heavy rain. Nor must they be dried in microwave ovens or other ovens.

Wipe the batteries carefully if moisture is present as it may affect their performance. Using an anti-humidity kit can help to avoid these problems and may even extend the life of your hearing instrument. Consult your Hearing Care Professional for further advice.

The chemicals in cosmetics, hairspray, perfume, after shave lotion, suntan lotion and insect repellent can damage your hearing instrument. You should always remove your hearing instrument before applying such products and allow time for the product to dry before reinserting your hearing instrument. If you use lotion, be sure to wipe your hands dry before putting on your hearing instrument.

## Using your hearing instrument

It takes time to adjust to a new hearing instrument. How long this adjustment takes differs from person to person. It will depend on a number of factors, such as whether you have had a hearing instrument before and the degree of your hearing loss.

### Seven easy steps to better hearing

#### 1. In the quiet of your home

Try to accustom yourself to all the new sounds. Listen to the many background sounds and try to identify each sound. Bear in mind that some sounds will seem different from what you are used to. You may have to learn to identify them again. Note that in time you will get accustomed to the sounds in your environment - if not please contact your Hearing Care Professional.

If using the hearing instrument makes you tired, take them off for a little while and have a rest. Gradually, you will begin to be able to listen for longer periods of time. Soon, you will be able to wear your hearing instrument comfortably all day long.



## 2. Conversation with another person

Sit with someone else in a quiet room. Face each other so you can read facial expressions easily. You may experience new speech sounds, which can seem a little disturbing in the beginning. However, after the brain has adapted to the new speech sounds, you should hear speech clearer.

## 3. Listen to radio or TV

When listening to the TV or the radio, start out by listening to news commentators since they usually speak clearly, then try other TV programs.

If you find it difficult to listen to TV or radio, ask your Hearing Care Professional for more information regarding the Oticon ConnectLine system\* and other available accessories.

## 4. In group conversations


Group situations are usually accompanied by a greater degree of background noise, and are, therefore, naturally more difficult to cope with. In such situations, therefore, focus your attention on the person you want to hear. If you miss a word, ask the speaker to repeat.

## 5. Telecoil use in church, theatre, or cinema

An increasing number of churches, theatres and public buildings often have loop systems installed. These systems send out wireless sound to be received by the telecoil in your hearing instrument. Typically, a sign will let you know whether the place has a loop system. Ask your Hearing Care Professional for details.

## 6. Using the telephone

When using the telephone tilt the receiver edge lightly on your cheek-bone and position the phone slightly over the ear close to the hearing instrument's microphone. The sound then flows directly into the hearing instrument's microphone opening. This way, the hearing instrument will not whistle and you ensure the best conditions to understand the conversation. When you have the receiver in this position, remember to speak directly into the microphone on the telephone in order to ensure good understanding in the opposite "end of the line".



If you find it difficult to use the telephone, ask your Hearing Care Professional for more information regarding the Oticon ConnectLine system\* and other available accessories.

If your hearing instrument has a telecoil (and your telephone has a built in tele-loop) you can switch into telecoil program in order to improve the sound reception further.

Remember that the telecoil in hearing instrument may pick up disturbing signals from electronic devices, such as a fax machines, computers, televisions or similar. Make sure that the hearing instrument has a distance of 2-3 meters to such devices when using the telecoil program.

### **7. Wireless and mobile phones**

Your hearing instrument is designed to comply with the most stringent Standards of International Electromagnetic Compatibility. However not all mobile phones are hearing instrument compatible. The varying degree of disturbance can be due to the nature of your particular mobile phone.

If you find it difficult to obtain a good result while using your mobile telephone, ask your Hearing Care Professional for more information regarding the Oticon ConnectLine system\* and other available accessories.

### **Use your hearing instruments all day long**

The best way to ensure better hearing is to practice listening until you can wear your hearing instruments comfortably all day. In most cases, infrequent use of a hearing instrument does not give you the full benefit. In most cases, you will not get the full benefit of the hearing instruments if you use them infrequently.

Your hearing instruments will not restore normal hearing. Nor will they prevent or improve a hearing impairment resulting from a physiological condition. What they will give you, however, is help towards making better use of the hearing ability that you have. If you have two hearing instruments always wear both.

*\*Not available for CIC, MIC and ITC Power instruments.*

<b>Common problems and their solutions</b>			
<b>Symptom</b>	<b>Possible causes</b>	<b>Solutions</b>	
<b>No sound</b>	Worn-out battery	Change battery	pg. 15
	Clogged sound outlet	Clean sound outlet or change ProWax	pg. 34 & 36
<b>Intermittent or reduced sound</b>	Clogged sound outlet	Clean sound outlet or change ProWax	pg. 34 & 36
	Moisture	Wipe battery and instrument with dry cloth	pg. 40
	Worn-out battery	Change battery	pg. 15
<b>Squealing noise</b>	Hearing instrument inserted improperly	Reinsert hearing instrument	pg. 19
	Ear wax accumulated in ear canal	Have ear canal examined by your doctor	

*If none of the above solutions solves the problem, ask your Hearing Care Professional for assistance.*

## Wireless accessories

### ConnectLine

As an enhancement to your wireless hearing instrument(s) a range of devices can be connected. The ConnectLine products can connect your hearing instrument(s) to a number of media and sound sources. This will allow you to receive audio signals from TV, phones, music players, PC or an external microphone wireless through your hearing instrument(s).

For detailed information on the ConnectLine system ask your Hearing Care Professional or visit [www.oticon.com](http://www.oticon.com)

### Not available for CIC, MIC and ITC Power instruments

## International Warranty

Your hearing instruments are covered by a limited warranty issued by the manufacturer for a period of 12 months from the date of delivery. This limited warranty covers manufacturing and material defects in the hearing instrument itself, but not accessories such as batteries, tubing, earwax filters etc.

Problems arising from improper handling or care, excessive use, accidents, repairs made by an unauthorized party, exposure to corrosive conditions, physical changes in your ear, damage due to foreign objects entering the device, or incorrect adjustments are NOT covered by the limited warranty and may void it.

The above warranty does not affect any legal rights that you might have under applicable national legislation governing sale of consumer goods. Your Hearing Care Professional may have issued a warranty that goes beyond the clauses of this limited warranty. Please consult him/her for further information.

### If you need service

Take your hearing instrument to your Hearing Care Professional, who may well be able to sort out minor problems and adjustments on the spot.

## Warranty Certificate

Name of Owner: \_\_\_\_\_

Dispenser: \_\_\_\_\_

Dispenser Address: \_\_\_\_\_

Dispenser Phone: \_\_\_\_\_

Purchase Date: \_\_\_\_\_

Warranty Period: \_\_\_\_\_ Month: \_\_\_\_\_

Model Left: \_\_\_\_\_ Serial no.: \_\_\_\_\_

Model Right: \_\_\_\_\_ Serial no.: \_\_\_\_\_

Battery Size: \_\_\_\_\_

## Mobile Phone

“Some hearing aid users have reported a buzzing sound in their hearing aids when they are using mobile phones, indicating that the mobile phone and hearing aid may not be compatible. According to the ANSI C63.19 standard (ANSI C63.19-2006 American National Standard Methods of Measurement of Compatibility Between Wireless Communications Devices and Hearing Aids), the compatibility of a particular hearing aid and mobile phone can be predicted by adding the rating for the hearing aid immunity to the rating for the mobile phone emissions. For example, the sum of a hearing aid rating of 2 (M2/T2) and a telephone rating of 3 (M3/T3) would result in a combined rating of 5. Any combined rating that equals at least 5 would provide “normal use;” a combined rating of 6 or greater would indicate “excellent performance.”

	<b>M</b>	<b>T</b>
<b>CIC, MIC</b>	M4	-
<b>CIC P, MIC P</b>	M4	-
<b>ITE, ITC</b>	M4	T3
<b>ITC P</b>	M2	T2

### **IMPORTANT NOTICE**

The performance of individual hearing instruments may vary with individual mobile phones. Therefore, please try this hearing instrument with your mobile phone or, if you are purchasing a new phone, be sure to try it with your hearing instrument prior to purchase. For additional guidance, please ask your mobile phone provider for the booklet entitled "Hearing Aid Compatibility with Digital Wireless Mobile Phones."

## **Technical information**

The hearing instrument contains a radio transmitter (*not applicable to CIC, MIC and ITC Power instruments*) using short range magnetic induction technology working at 3.84 MHz. The magnetic field strength of the transmitter is  $< -42 \text{ dB}\mu\text{A/m @ 10m}$ .

The emission power from the radio system is well below international emission limits for Human Exposure. For comparison, the radiation of the hearing instrument is lower than unintended electromagnetic radiation from for example halogen lamps, computer monitors, dishwashers etc. The hearing instrument complies with international standards concerning Electromagnetic Compability.

Due to the limited size available on the instruments all relevant approval markings are found in this document.

The ITC/ITE hearing instruments contains a module with:

FCC ID: U28FUITE03

IC: 1350B-FUITE03

The device complies with Part 15 of the FCC rules and RSS-210 of Industry Canada.

Operation is subject to the following two conditions:

1. this device may not cause harmful interference.
2. this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Hereby, Oticon, declares that this hearing instrument is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC. Declaration of conformity is available at:

Oticon A/S  
Kongebakken 9  
DK-2765 Smørum  
Denmark  
[www.oticon.com](http://www.oticon.com)

CE 0543 0682



Waste from electronic equipment must be handled according to local regulations.



NI 175

## Settings overview for your hearing instrument

Hearing Instrument				
Left			Right	
<input type="checkbox"/> Yes	<input type="checkbox"/> No	Volume control	<input type="checkbox"/> Yes	<input type="checkbox"/> No
<input type="checkbox"/> Yes	<input type="checkbox"/> No	Program shift	<input type="checkbox"/> Yes	<input type="checkbox"/> No
<input type="checkbox"/> Yes	<input type="checkbox"/> No	Mute	<input type="checkbox"/> Yes	<input type="checkbox"/> No
Volume control indicators				
<input type="checkbox"/> On	<input type="checkbox"/> Off	Beeps at min /max volume	<input type="checkbox"/> On	<input type="checkbox"/> Off
<input type="checkbox"/> On	<input type="checkbox"/> Off	Clicks when changing volume	<input type="checkbox"/> On	<input type="checkbox"/> Off
<input type="checkbox"/> On	<input type="checkbox"/> Off	Beeps at preferred volume	<input type="checkbox"/> On	<input type="checkbox"/> Off
Battery indicators				
<input type="checkbox"/> On	<input type="checkbox"/> Off	Low battery warning	<input type="checkbox"/> On	<input type="checkbox"/> Off

Hearing Program Description	
Program No.	
1	
2	
3	
4	



908 73 410 00 / 04.12

908 73 410 00 / 04.12

## People First

People First is our promise  
to empower people  
to communicate freely,  
interact naturally and  
participate actively



9087341000000001

**oticon**  
PEOPLE FIRST