

# Diatec Consumables Catalogue 2020



# Contents



Page 4 – 7

## Disposable ear tips



Page 8 – 9

## Probe tips



Page 10 – 12

## Supplies for insert earphones



Page 13 – 14

## ABR & newborn



Page 15 – 16

## Hearing Aid Fitting & supplies



Page 17 – 18

## Hearing Aid Maintenance & hearing care



Page 19 – 21

## Impression materials & supplies



Page 22 – 26

## Disposable electrodes



Page 27 – 28

## Electrode cables



Page 29 – 30

## Patient preparation



Page 31 – 34

## Infection control



Page 35

## Balance supplies



Page 36 – 37

## Otoscopy Supplies



Page 38

## Clinical tools



Page 39 – 41

## Paper



Page 42 – 44

## Wax removal accessories



Page 45 – 49

## Parts for audiometry



Page 50

## Miscellaneous

### List of abbreviations

BAL

Balance

HIF

Hearing Instrument Fitting

IMP

Impedance

ABR

Auditory Brainstem Response

AUD

Audiometry

OAE

Otoacoustic Emissions



Single patient use only!



Product items are not made with natural rubber latex

# Selecting the correct ear tip

Age	Ear tips inserted against ear canal	Ear tips inserted against entrance of ear canal			Ear tips for Echo-Screen®	Ear tips for AccuScreen®			
0.5 y/o			Green pediatric		 4.0 mm, ANA	 3.7 mm, AOT			
			3-5 mm, flanged		3-6 mm, green		6 mm, green	 4.5 mm, ANA	 4.0 mm, AOT
	 7.5 mm, umbrella		4-7 mm, flanged		6 mm, yellow	 infant tree, ANA	 5.0 mm, AOT		
	 9 mm, umbrella		5-8 mm, flanged			 mushroom, ANA	 4-7 mm, AOT		
	 10 mm, umbrella		7 mm, mushroom		7 mm, red		7 mm, yellow		
	 11 mm, umbrella		8 mm, mushroom		8 mm, blue		8 mm, red		
	 13 mm, umbrella		9 mm, mushroom		9 mm, green		9 mm, blue		
			10 mm, mushroom		10 mm, yellow		10 mm, green		
	 15 mm, umbrella		11 mm, mushroom		11 mm, red		11 mm, yellow		
			12 mm, mushroom		12 mm, blue		12 mm, red		
Adults	 19 mm, umbrella		13 mm, mushroom		13 mm, green		13 mm, blue		
			14 mm, mushroom		14 mm, yellow		14 mm, green		
			15 mm, mushroom		16 mm, red		16 mm, yellow		
	 22 mm, umbrella		19 mm, mushroom		18.5 mm, blue				



Selecting the correct size of ear tip takes practice

In general, flanged or mushroom-shaped ear tips will be used for most testing (Impedance, DPOAE, TEOAE and Automated ABR).

Sometimes, an umbrella-shaped ear tip will be used for performing a quick impedance screening.

The table can be used as a guide to help you select the appropriate ear tip size.

## Disposable Sanibel™ ADI ear tips

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip: Amplivox: all OtoWave versions. Interacoustics: Titan, AT235, AA220, Eclipse, MT10, Sera. Maico: EasyTym, MI 22, MI 23, MI 24, MI 26, MI 34, MI 44, Easyscreen



NEW

### ADI Assortment Box

Part no.	Description	Unit
8511445	Various sizes for Tympanometry and OAE	
8501572	Screening ear tip kit, various sizes for newborn screening	



### Flanged infant ear tips

Part no.	Description	Unit
8502838	Pediatric (green)	Bag/100 pcs
8012964	3-5 mm (red)	Bag/100 pcs
8012966	4-7 mm (blue)	Bag/100 pcs
8012968	5-8 mm (yellow)	Bag/100 pcs



### Umbrella ear tips

Part no.	Description	Unit
8012986	7.5 mm (green)	Bag/100 pcs
8012988	9 mm (red)	Bag/100 pcs
8012990	10 mm (blue)	Bag/100 pcs
8012992	11 mm (yellow)	Bag/100 pcs
8012994	13 mm (green)	Bag/100 pcs
8012996	15 mm (red)	Bag/100 pcs
8012998	19 mm (blue)	Bag/100 pcs
8013000	22 mm (yellow)	Bag/100 pcs



### IA series ear tips

For use with many acoustic immittance instruments and otoacoustic emission systems, including: GSI: TymStar Pro, Allegro, Novus and Corti. Interacoustics (clear probe): AT235, AA222 and Titan and Otoread II, MT10. Maico: Touchtym Racecar, MI24, MI26, MI34, MI36, Easytym. Madsen/Otometrics: Zodiac (new), Accuscreen. Path Medical: all Sentiero products,



### Infant ear tips

Part no.	Description	Unit
8013006	3 mm (yellow)	Bag/100 pcs
8002035	4 mm flanged (red)	Bag/100 pcs



### Mushroom ear tips

Part no.	Description	Unit
8013002	7 mm (blue)	Bag/100 pcs
8013004	8 mm (yellow)	Bag/100 pcs
8012970	9 mm (green)	Bag/100 pcs
8012972	10 mm (red)	Bag/100 pcs
8012974	11 mm (blue)	Bag/100 pcs
8012976	12 mm (yellow)	Bag/100 pcs
8012978	13 mm (green)	Bag/100 pcs
8012980	14 mm (red)	Bag/100 pcs
8012982	15 mm (blue)	Bag/100 pcs
8012984	19 mm (yellow)	Bag/100 pcs



### Mushroom ear tips

Part no.	Description	Unit
8520348	3-6 mm (green)	Bag/100 pcs
8520349	6 mm (yellow)	Bag/100 pcs
8520350	7 mm (red)	Bag/100 pcs
8520351	8 mm (blue)	Bag/100 pcs
8520352	9 mm (green)	Bag/100 pcs
8520353	10 mm (yellow)	Bag/100 pcs
8520354	11 mm (red)	Bag/100 pcs
8520355	12 mm (blue)	Bag/100 pcs
8520356	13 mm (green)	Bag/100 pcs
8520357	14 mm (yellow)	Bag/100 pcs
8520358	16 mm (red)	Bag/100 pcs
8520359	18.5 mm (blue)	Bag/100 pcs



Using a new ear tip is optimal for obtaining a seal when inserting it into the patients' ear. A multi-used ear tip may cause inaccuracies during testing.



# Disposable ear tips

## AKR Mushroom ear tips

Use with the following acoustic immittance instruments and otoacoustic emission systems with BLACK probe tip from Interacoustics and Maico: AT235, AA220, AA220, Eclipse, MT10, Otoread (I), RaceCar MI 22, MI 23, MI 24, MI 26, MI 34, MI 44. Madsen/Otometrics: Capella 2, Zodiac 901, Biologic AuDX Pro (+Flex). Neurosoft: Neuro-Audio, Oscilla: TSM 300/400 + T800 series. All GSI tympanometers from GSI 20 - 39 + Tymptstar.



### Mushroom ear tips

Part no.	Description	Unit
8518409	AKR Assortment Box, various sizes	box
8500570	3-6 mm (yellow)	bag/100 pcs
8107510	4 mm (red)	bag/100 pcs
8107511	5 mm (blue)	bag/100 pcs
8500560	6 mm (green)	bag/100 pcs
8500550	7 mm (yellow)	bag/100 pcs
8500540	8 mm (red)	bag/100 pcs
8500530	9 mm (blue)	bag/100 pcs
8500520	10 mm (green)	bag/100 pcs
8500510	11 mm (yellow)	bag/100 pcs
8500500	12 mm (red)	bag/100 pcs
8500490	13 mm (blue)	bag/100 pcs
8107520	14 mm (green)	bag/100 pcs
8500480	16 mm (yellow)	bag/100 pcs

NEW



ABR

OAE

## DISPOSABLE SANIBEL™ AOT and ANA ear tips



### ANA ear tips for echo-Screen®

For use with Natus: echo-Screen®, Madsen/Otometrics: AccuScreen® Classic. Interacoustics: Eclipse with pink probe tip

Part no.	Description	Unit
8106363	4.0 mm (blue)	Bag/50 pcs
8106362	4.5 mm (red)	Bag/50 pcs
8106364	5.0 mm (green)	Bag/50 pcs
8106366	Infant tree (yellow)	Bag/50 pcs
8106365	Mushroom (green)	Bag/50 pcs



ABR

OAE

## IMPEDANCE EAR TIPS – Fading product

For use with the following systems with BLACK probe tip: Old MT10, AT235, AA22x, MI 22, MI 23, MI 24, MI 26, MI 34, MI 44. FADING OUT product: consider AKR instead.



### Umbrella ear tips

Part no.	Description	Unit
8107617	10 mm (green)	bag/100 pcs
8107618	13 mm (green)	bag/100 pcs
8107619	16 mm (green)	bag/100 pcs
8107620	19 mm (blue)	bag/100 pcs
8107621	22 mm (blue)	bag/100 pcs



IMP



### Sanibel™ AZE mushroom ear tips

For use with old MT10 and the following systems with black probe tip: AT235, AA220, MI22, MI23, MI24, MI26, MI34, MI44, and old SPL60 (4mm)

Part no.	Description	Unit
8107622	6 mm (violet)	Bag/100 pcs
8107623	7 mm (yellow)	Bag/100 pcs
8107625	10 mm (blue)	Bag/100 pcs
8107626	11 mm (white)	Bag/100 pcs
8107629	18 mm (blue)	Bag/100 pcs



IMP



### AOT ear tips

For use with Madsen/Otometrics AccuScreen®, Biologic Bio-logic® AuDX®, Pathmedical: all Sentiero products

Part no.	Description	Unit
8106367	3.7 mm (yellow)	Bag/50 pcs
8106368	4.0 mm (light green)	Bag/50 pcs
8106369	4.5 mm (light pink)	Bag/50 pcs
8106370	5.0 mm (light blue)	Bag/50 pcs
8106371	4-7 mm (transparent)	Bag/50 pcs



ABR

OAE

## Grason & Associates ot & MO select series ear tips



### OT select ear tips w/ triangle holes

Ear tips with triangle holes for use with Otodynamics UGD, SND, SGD probes

Part no.	Description	Unit
8102644	4.8 mm (blue)	Bag/100 pcs
8102645	6 mm (green)	Bag/100 pcs
8102646	8 mm (red)	Bag/100 pcs
8102647	10 mm (blue)	Bag/100 pcs
8102648	12 mm (green)	Bag/100 pcs
8102649	14 mm (red)	Bag/100 pcs
8102650	16 mm (blue)	Bag/100 pcs



ABR

OAE



### OT select ear tips w/two probe holes

Ear tips with two probe holes for use with Otodynamics ECP, UGS, SNS, SGS probes.

Part no.	Description	Unit
8102632	3.5 mm (blue)	Bag/100 pcs
8102633	4.5 mm (blue)	Bag/100 pcs
8102634	5.5 mm (blue)	Bag/100 pcs
8102635	6.5 mm (blue)	Bag/100 pcs
8102636	4.8 mm (blue)	Bag/100 pcs
8102637	6 mm (green)	Bag/100 pcs
8102638	8 mm (red)	Bag/100 pcs
8102639	10 mm (blue)	Bag/100 pcs
8102640	12 mm (green)	Bag/100 pcs
8102641	14 mm (red)	Bag/100 pcs
8102642	16 mm (blue)	Bag/100 pcs



ABR

OAE



### MO select ear tips for Madsen OTOflex

For MADSEN® OTOflex™.

Part no.	Description	Unit
8102653	3-6 mm (yellow)	Bag/100 pcs
8102654	6 mm (red)	Bag/100 pcs
8102655	7 mm (blue)	Bag/100 pcs
8102656	8 mm (green)	Bag/100 pcs
8102657	9 mm (yellow)	Bag/100 pcs
8102658	10 mm (red)	Bag/100 pcs
8102659	11 mm (blue)	Bag/100 pcs
8102660	12 mm (green)	Bag/100 pcs
8102661	13 mm (yellow)	Bag/100 pcs
8102662	14 mm (red)	Bag/100 pcs
8102663	16 mm (blue)	Bag/100 pcs
8102664	18.5 mm (green)	Bag/100 pcs



ABR

OAE



### Otodynamics Ear Tips

For UGD or SND probes

Part no.	Description	Unit
T-R4.8S	4.8mm (blue)	Bag/100 pcs
T-R5.8B	5.8mm (blue)	Bag/100 pcs
T-R6.8B	6.8 mm (blue)	Bag/100 pcs



### Otodynamics Ear Tips

Eartips for measuring DPOAE and TEOAE. For use with Otodynamics ILO-OAE instruments with DPOAE and TEOAE probes.

Part no.	Description	Unit
T-R7M	7mm (blue)	Bag/100 pcs
T-R8M	8mm (blue)	Bag/100 pcs
T-R9M	9mm (blue)	Bag/100 pcs
T-R11M	11mm (blue)	Bag/100 pcs
T-R13M	13mm (blue)	Bag/100 pcs

### Ear tips for OAE screeners

For use with Biologic AuDx, Welch Allyn OAE screener and Natus echo-Screen® III.



### OAE ear tips, AuDx

Part no.	Description	Unit
8107885	AuDx Premie 3 mm	Bag/100 pcs
8107889	AuDx Tree 3.5-5mm	Bag/100 pcs
8107887	AuDx Infant 6 mm	Bag/100 pcs
8102681	Adult OAE foam	Bag/100 pcs
8102682	Pediatric OAE foam	Bag/100 pcs
8103123	Jumbo OAE foam	Bag/100 pcs



ABR

OAE





# Disposable ear tips



## Grason-Stadler eartips

Compatible with Tympan, GSI33, GSI38, and GSI39 when testing Contra.

Part no.	Description	Unit
8012929	7mm (yellow)	Box/12 pcs
8012930	8mm (pink)	Box/12 pcs
8012931	9mm (blue)	Box/12 pcs
8012932	10mm (green)	Box/12 pcs
8012933	11mm (pink)	Box/12 pcs
8012934	12mm (yellow)	Box/12 pcs
8012935	13mm (blue)	Box/12 pcs
8012936	14mm (green)	Box/12 pcs



LATEX

OAE

## Grason & Associates mixed ear tips



### TS211

Smaller offices needing the full size selection of 4-16mm may prefer the TS211 with 10 of each KR ear tip, sizes 4-16mm.



### TS291

Smaller offices needing the full size selection of 3-14mm may prefer the TS291 with 10 of each GSI mushroom ear tip, sizes 3-14mm.



## Clinical Audiology Kit

Suitable for students, graduates and new clinics

Part no.	Description	Unit
1100141		1 kit

Includes:  
 1x Heine Mini 3000 Otoscope  
 1x LED Earlight  
 50x Heine Disposable Specula 4.0mm  
 50x Heine Disposable Specula 2.5mm  
 1x Stetaset  
 1x Combo 3 in 1 tool  
 1x Vent cleaner  
 5x Sanibel Universal Type A Probe tube  
 20x Alcohol Pads  
 20x 10mm Green impedance tips  
 20x 11mm Yellow impedance tips  
 20x 12mm Red impedance tips  
 20x 13mm Blue impedance tips  
 20x 14mm Green impedance tips  
 20x 3A foam tips  
 20x 3B foam tips  
 20x Disposable ear cushion covers  
 2x Clear plastic containers



## DP Probe ear tips

Large foam tips for adult ears.

Part no.	Unit
T-R100L	Bag/100 pcs





## Probe tips



IMP

### Probe tip for Sanibel™ ADI™ ear tips

Fits Otowave and new MT10™.

Part no.	Unit
8002592	1 pc



ABR IMP OAE

### Probe tip for Sanibel™ ADI™ ear tips

Fits Titan™ and easyTymp™.

Part no.	Unit
8001701	1 pc




IMP

### Probe tip kit

Fits AT235™ or newer, AA222™, Sera™, Tympstar Pro™ Novus™, easyScreen™, and touchTymp™. Includes probe tip and gasket.

Part no.	Unit
8107762	Set/2 pcs




100 pcs

ABR OAE

### Sanibel™ probe tubes for ADI™ ear tips

Fits Corti™, ERO•SCAN™ and the new OtoRead™.

Part no.	Unit
8104159	Box/100 pcs




ABR OAE

### Probe Tube Removal Tool

Fits Corti™, ERO•SCAN™ and the new OtoRead™.

Part no.	Unit
8513131	Bag/1 pc



+

1 pc +  
100 pcs

ABR OAE

### Probe Tube Removal Tool kit

Fits Corti™, ERO•SCAN™ and the new OtoRead™.

Part no.	Unit
8513132	Bag/100 pcs + 1 pc



OAE

### Probe tip kit for OAE Ø3.3 mm ear tips

Fits all Eclipse™ – use with ADI™ ear tips. Requires no re-calibration.

Part no.	Unit
8012960	Bag/8 pcs



IMP

### Sanibel™ probe tip adaptor kit

Converts AT235™, AT235h™, AA220™, AA222™(blue/new), MI24™, MI26™, MI34, MI44™ for use with ADI™ ear tips. Includes probe tip, washer, damper. No recalibration required.

Part no.	Unit
8500430	Set/3 pcs



OAE

### Remote probe tip replacements for ERO•SCAN™

Fits ERO•SCAN™.

Part no.	Unit
8102340	Bag/4 pcs



## Probe tips



OAE

### Internal probe tip replacements for ERO•SCAN™

Fits ERO•SCAN™.

Part no.	Unit
8102339	Bag/4 pcs



OAE

### ER10D probe tips for ERO•SCAN™ Pro and OtoRead™

It is recommended to clean the probe tip after every patient and to change the probe tip after 20 tests (40 ears). No calibration necessary.

Part no.	Unit
8012836	Bag/4 pcs



OAE

### Internal probe tip replacements for ERO•SCAN™ Pro

Fits ERO•Scan™ Pro.

Part no.	Unit
8110364	Bag/4 pcs

## Sanibel™ Probe Tube Removal Tool

A = Probe

B = Probe Tube Installer

C = Attachment to probe cable


D = Tube Remover

To Remove a Probe Tube

1. Place the front hole of the Probe Tube Tool over the end of the probe, as shown. The probe should be seated against the face of the tool and snap in place.
2. Squeeze the tool closed and hold it.
3. Twist the tool a couple of times while holding the tool closed and pulling away from the probe.
4. The probe tube will pull out from the probe. Discard the probe tube.

To Install a New Probe Tube

5. Place a new tube into the hole on the top of the tool.
6. Align the probe with the tube.
7. Push the probe onto the probe tube until the probe contacts the face of the probe tube tool and a snapping sound is heard.
8. Slide the tool off leaving the tube in place in the probe.





**3A E-A-RLINK™ foam ear tips**  
Complies to all insert phones. 13.7mm

Part no.	Description	Unit
8106226	Adult/Medium size	Bag/100 pcs



**3B E-A-RLINK™ foam ear tips**  
Complies to all insert phones. 9.7mm

Part no.	Description	Unit
8106227	Pediatric/Small size	Bag/100 pcs



**3C E-A-RLINK™ foam ear tips**  
Complies to all insert phones.

Part no.	Description	Unit
8500100	Jumbo/Large size 17.8 mm	Bag/24 pcs




**E-A-RLINK™ infant ear tips**  
For use with infants or in very small ear canals.

Part no.	Description	Unit
8500300	3.5 mm	Bag/20 pcs
8500310	4.0 mm	Bag/20 pcs




**Sanibel™ ADI ear tip adaptor**  
Makes Sanibel™ ADI ear tips compatible with 3A insert phones.

Part no.	Unit
8500580	Bag/6 pcs




**Sound tubes for E-A-Rtone® 3A ABR 50 ohm**  
Includes adaptor nipples, 1 pc blue and 1 pc red 50 ohm.

Part no.	Unit
8500720	Set/2 pcs



**Sanibel™ IP30™ tube connectors**  
For use with IP30™ and E-A-Rtone® 3A.

Part no.	Unit
8107391	Bag/10 pcs




**Tube nipples for E-A-Rtone® 3A**  
For use with sound tubes 3A.

Part no.	Unit
8500140	Bag/10 pcs



**Sanibel™ sound tubes and Sanibel™ ADI™ ear tip adaptor kit**  
Fits E-A-Rtone® 3A inserts and IP30™ insert phones. Includes ADI™ ear tip adaptors, 1 pc blue and 1 pc red 50 ohm.

Part no.	Unit
8100590	Set/2 pcs



### Sound tubes for E-A-Rtone® 3A 10 ohm

Includes adaptor nipples.

Part no.	Unit
8107396	Bag/4 pcs



### Sound tubes for E-A-Rtone® 3A 10 ohm

Includes adaptor nipples.

Part no.	Unit
8103605	Bag/4 pcs



### Sanibel™ IP30 and E-A-Rtone® sound tubes and nipples

Fit all IP30™ and E-A-Rtone® versions. Includes tube connectors, 1 pc blue and 1 pc red 50 ohm.

Part no.	Unit
8501072	Set/2 pcs



### E-A-Rtone® replacement fastener clip

Clip for fastening insert phone on clothes.

Part no.	Unit
8000744	1 pc



### Sanibel™ fastener kit for IP30™

Includes velcro neckstrap and clip.

Part no.	Unit
8507410	Set



### Sanibel™ ear tip adaptor cleaning tool

Small brush for cleaning the Sanibel™ probe tip adaptors before disinfection.

Part no.	Unit
8500830	Bag/10 pcs



ABR

### ABR consumable kit

Includes electrodes, insert phone tips, skin scrub and gazue.

Part no.	Description	Unit
8109670	Adult 3A	Kit/100 tests
8109671	Pediatric 3B	Kit/100 tests



HIF

### Wax Brush 3 in 1

Includes:

- Full brush to remove unwanted wax from hearing aids.
- Magnet at the top of the brush to remove batteries
- Wire look at opposite end

Part no.	Unit
9003569	1 pc



HIF

### Temporary Moulds

Used as a temporary/test tip with any Post Auricular Hearing Aid (BTE). Pre-bent Clear #13 tube. Standard 13.7mm diameter.

Part no.	Unit
9002602	Bag/25 pcs




ABR

### Sanibel™ Infant EarCup™ adaptor set

Connects insert phones with Infant EarCups™.  
NB: Do not contain filters.

Part no.	Unit
8500450	Set/2 pcs



ABR

### Sanibel™ sound tubes and Sanibel™ Infant EarCup™ adaptor kit

Includes adaptors, 1 blue and 1 red tube.

Part no.	Unit
8100577	Set/2 pcs



### 3M Coban Self-Adherent Wrap

100mm x 2m. Sticks to itself without need for adhesive, pins or clips for fast and easy application.

Part no.	Unit
9000420	1 Roll



ABR

### Newborn wipes

Fragrance free.

Part no.	Unit
9000287	Canister/ 125 wipes



### (DPC) DP Probe Couplers

For use with UGD, SND, SGD probes

Part no.	Unit
DPC-10	Bag/10 pcs
DPC-100	Bag/100 pcs





2

ABR

## Sanibel™ Infant EarCup™ with snap or tab electrodes

Reduce test time! Test on both ears in one screening. Compatible with Titan™, Eclipse™, AccuScreen® (1st generation), MB II Classic and Echo-Screen®. Use only with equipment calibrated for over-ear style disposables. Disposable over-the-ear earphones are high quality, low cost, and ideal for bilateral newborn hearing screenings. 1 set consists of 2 Infant EarCups™, 3 snap or tab electrodes and 1 alcohol pad.

For use with a range of different providers of newborn hearing screening equipment:

Grason Stadler: Novus. Interacoustics: Titan, Sera. Maico: Easyscreen, MBII Classsic. Madsen/Otometrics: Accuscreen. Natus: Biologic Abear, Echoscreen.

Part no.	Description	Unit	
8049600	Tab	Box/20 sets	
8049700	Snap	Box/20 sets	
8508848	w/o electrodes	Box/20 sets	\$245

## Sanibel Infant EarCup™ w/ Hydrogel.

New Hydrogel version to ensure a more comfortable solution for the baby while testing.


Choose any electrode type you wish to add to meet your preference.

Order electrodes separately

Part no.	Unit
8516095	Box/20 sets

NEW





The "REM effect"  
– better fitting & fewer follow-ups!

## Use Sanibel™ probe tubes for any equipment

We often hear from the industry that only the manufacturers' specific probe tubes should be used with their specific equipment. However, a recent comparative study by Sanibel Supply® disproves this myth.

In this study, our in-house audiologists and engineers compared the effects of different probe tubes on an Affinity<sup>2.0</sup> system. They compared AudioScan tubes (red), Otometrics tubes (blue), Sanibel™ tubes (orange) and Sanibel™ tubes (grey).

This involved running each tube via the 'tube calibration' process and then comparing the following REUG curve when the tube had been left in the calibration position.

The investigation showed no significant difference in the REUG curve and all probe tubes achieved the same result when the calibration process had been run correctly.



50 pcs

HIF

### Silicone tubes for SPL60™

For use with Affinity™ REM probe.

Part no.	Unit
8500700	Bag/50 pcs




HIF

### Sanibel™ REM probe tubes, type A

Universal fit, length 76 mm, beveled 45° angle tip and black marker ring.

Part no.	Unit
8104116	Tray/50 pcs



HIF

### Sanibel™ REM probe tubes, type B

Cost efficient, length 80 mm and marker ring.

Part no.	Unit
8500360	Bag/40 pcs



HIF

### Audioscann probe tubes

Compatible with Audioscan, Axiom and Verifit equipment

Part no.	Unit
RE407-40	Tray0 pcs





HIF

## Pre-bent tubing

Part no.	Description	Unit
8515641	Standard #13	Bag/25 pcs
8515642	Standard #12	Bag/25 pcs
8515643	Thick #13	Bag/25 pcs
8515644	Extra Thick #13	Bag/25 pcs

2  
HIF

## Tube lock with tubing

Brass metal lock

Part no.	Description	Unit
8515658	Thick	Bag/10 pcs
8515659	Super Thick	Bag/10 pcs



## Surgical Blades and Handles

Suitable for cutting HA tubing.

Part no.	Description	Unit
8515427	Surgical blade handle No.3	1 pc
8515428	Surgical blades No11	100 pcs




New

## Tubing Puller

For replacing worn out tubing, Pulls tubes easily through any type of mold or canal tip

Part no.	Unit
9004935	1 pc



HIF

## Pre-bent tubing

Tinted

Part no.	Description	Unit
8515652	Standard #13	Bag/10 pcs
8515653	Medium #13	Bag/25 pcs
8515644	Thick #13	Bag/25 pcs

2  
HIF

## Tube Lock without tubing

Brass Tube-Locks Tube-Locks are used for attaching tubing to flexible type earmolds without requiring any cement, yet allowing tubing replacement as needed.

Part no.	Description	Unit
8515664	#13 Thick	1 pc
8515665	#13 Super Thick	1 pc



New

## Tube lock applicator

A handy tool when inserting Tube-Lock into mold.

Part no.	Unit
9002907	1 pc

HIF

## Loctite 406 Instant adhesive

An instant adhesive designed especially for the fast bonding of plastics and elastomeric materials where very fast fixturing is required. This product offers exceptionally low viscosity.

Part no.	Unit
8005067	1 unit







## Take care of your hearing aids

Take good care of your hearing aid. Your Hearing instrument is a high tech device needing regular care in order to provide you the best performance every day.



New

### Cleansing tablets

Part no.	Unit
8520148	4 pk



New

### Dehydrating capsules

Part no.	Unit
8520149	20 pk



New

### Drying cup

Part no.	Unit
8520150	1 pc



New

### Cleaning cup

Part no.	Unit
8520151	1 pc



New

### Otofloss

Part no.	Unit
8520152	1 pc



New

### Perfect Dry LUX box for all types of hearing aids

Part no.	Unit
8520154	1 pc



New

### Starter kit

Part no.	Unit
8520155	kit



New

### Drying spray

Part no.	Unit
8520153	1 pc



New

### Cleaning spray

Part no.	Unit
8520147	1 pc



New

### Cleaning wipes

Part no.	Unit
8520146	30 pk



New

### Protect your ears

Part no.	Description
8520052	Natural lubricant



New

### Protect your ears

Part no.	Description
8520053	Dermoprotective Gel



New

### Protect your ears

Part no.	Description
8520063	Earspray



New

### Sanohra Earplugs Fly

Sanohra Fly ear plugs. Reduces ear pressure experiences when flying.

Part no.	Description	Unit
8519351	EarPlugs Fly Kids, EN	1 set
8519352	EarPlugs Fly Adult, EN	1 set



### Sanohra Earplugs Swim

Sanohra Swim earplugs. Enables people to swim without getting water in the ear, reducing ear infections.

Part no.	Description	Unit
8519353	EarPlugs Swim Kids, EN	1 set
8520143	EarPlugs Swim Adult, EN	1 set

### SANOHR A fly ear plugs

Beyond "noise protection" it can be important to protect the ear drum (tympanic membrane) from changes in external pressure that are too rapid, e.g., during take off and landing phases in airplanes. For airplane passengers experiencing a problem with the pressure equalisation, the use of a special ear protection product means freedom from an unpleasant feeling or piercing pain during flight.

### SANOHR A swim ear plugs

In case of injury to the ear drum or insertion of a vent tube, it is necessary to prevent water from reaching the ear drum through the outer auditory canal. Swimming, bathing or showering might cause inflammation of the middle ear. Staying in chlorinated water for an extended period of time can cause the skin in the outer auditory canal to swell and cause inflammation (otitis externa).





## Foam Otoblock™

Part no.	Description	Unit
Otblkpink	Small (pink)	Pkt/100 pcs
Otblkwhite	Medium (white)	Pkt/50 pcs
Otblklarge	Large (blue)	Pkt/50 pcs
825403202	Small (grey)	Pkt/50 pcs
8254003303	Medium (white)	Pkt/50 pcs
W05040	Large (blue)	Pkt/50 pcs



HIF

## Earmold Air Blower

Part no.	Unit
8515562	Box/12 pcs



## Cotton Otoblock™

Part no.	Description	Unit
9140797	Small (blue)	Pkt/50 pcs
9140798	Medium (white)	Pkt/50 pcs
9140799	Large (red)	Pkt/50 pcs



## Professional Syringes

Part no	Description	Unit
8900102709	4mm Grey	1 pc
923000002	3mm Blue	1 pc



HIF

## OTO-EASE™ Earmold & Hearing Aid Lubricant

A grease-less bacteria free liquid lubricant which eases insertion of ear-molds, I.T.E and C.I.C hearing aids

Part no.	Unit
9003259s	1x 15ml bottle



## Mini Cleaning Brush

Part no.	Unit
8251805605	1 pc



## Flat Type Spatula

Part no.	Unit
9000231	1 pc



## Bite Blocks

Made from styrofoam these disposable bite blocks are used to hold the patient's mouth open in a comfortable position when taking open jaw ear impressions

Part no.	Unit
8519356	Bag/100 pcs



## A-ZoFT Ear Impression Material

B & C components for ear impressions

Part no.	Unit
87825015105	40 foil bags
87825015004	2x 800g jars





New

## Dreve 4545 Otoform AK Impression Material Kit

8 x S 50 double cartridges of 2 x 25ml A+B

Part no.	Unit
9004439	2x25ml



New

## Otoform A/K 2 part silicone

Proven kneadable impression material, non-sticky and smooth with a slightly pressure-forming consistency. An even and controlled pressure is applied onto the ear tissue. Otoform AK is advantageous for little resilient ear tissues for fitting and sealing of earmolds

Part no.	Unit
9004438	2 pcs

New



## Mixing Pads

Disposable paper pads for mixing impression material. Foam, no slip surface on bottom of pad. Package of two pads with 35 sheets per pad. Dimensions 6" x 6"

Part no.	Unit
8522106	100 sheets

New




## Detax Pink Putty Impression Kit

This impression putty has a soft creamy consistency for easy mixing and removal. The Kit includes two dosing spoons and 500ml of catalyst and 500ml of base material to produce approximately 60 impressions.

Part no.	Unit
A16738	set

New



## Mixing tray

Comes with dosing scale

Part no.	Unit
8522105	1 pc



## Stethoclip

Hearing Instrument listening stethoscope for reasuring the fitting process

Part no.	Description	Unit
8515600	Stethoclip, black	1 unit

Large soft  
electrode size!

60 pcs



ABR

### Sanibel™ snap electrodes

Pre-gelled with a very gentle hydrogel, perfect for ABR testing on infants. Size 2.2 x 3 cm. Item no. 8107137 is packed in 20 smaller pouches of 3 electrodes.

Part no.	Unit
8107137	Box/60 pcs (20x3 pcs)



Large soft  
electrode size!

60 pcs



ABR

### Sanibel™ tab electrodes

Pre-gelled cloth electrode with a very gentle hydrogel, perfect for ABR testing on infants. Size 2.2 x 3 cm. Packed in 20 smaller pouches of 3 electrodes.

Part no.	Unit
8107410	Box/60 pcs (20x3 pcs)



60 pcs



ABR

### Sanibel™ wetgel snap electrodes

Pre-gelled round foam ECG electrode with button connector. Packaged in convenient 12x5 pcs re-sealable pouches.

Part no.	Unit
8108563	Bag/60 pcs (12x5 pcs)



100 pcs



ABR

### Sanibel™ vinyl tab electrode

For alligator clip cable. Solid gel. Size 2.5 x 2.2 cm. Packaged in convenient re-sealable pouches.

Part no.	Unit
8107426	Bag/100 pcs (10x10 pcs)



ABR

### Sanibel™ TM electrode

For ECoG. Use with cable 8500650.

Part no.	Unit
8011489	Box/10 pcs
8049550	1 pc



300 pcs

Small size and soft  
electrode!



ABR

### Sanibel™ prewired electrodes

Prewired electrodes with a small footprint. Very good for testing infants. Solid gel. Size 2.2 x 2.2 cm. 100 pouches of 3 pcs.

Part no.	Unit
8107425	Box/300 pcs



20 pcs



ABR

### E-A-RLINK™ gold foil tip electrodes

For ECoG. Use with cable 8011343.

Part no.	Description	Unit
8500370	Adult 13 mm	Bag/20 pcs
8500440	Pediatric 10 mm	Bag/20 pcs



NEW

50 pcs



ABR

### Vermed Economy line 3-pak Solid gel snap Electrode

Solid Gel type snap electrode suitable for newborn screening and Diagnostic follow-up. Latex free. Intended use areas: Newborn, Diagnostic follow-up all ages.

Part no.	Unit
8516252	Bag/50 pcs




ABR

### Ambu® Ground electrode 71410-M/1

Part no.	Unit
8500690	1 pc



# Disposable electrodes




ABR

## Vermed Wetgel Snap Electrode

Pre-gelled hydrogel electrode with Ag/AgCl sensor. Very suitable for newborn screening or diagnostic follow up. Packaged 3x20 pcs. per pouch.

Intended use areas: Newborn, Diagnostic follow-up, ABR, C-Vemp.

Part no.	Unit
8516221	Bag/60 pcs (3x20 pcs)




ABR

## Vermed Tab Electrode

Rectangular tab electrode for newborn screening. 12,7x22,9 mm. Conductive Adhesive Hydrogel. 10x10 pcs. per pouch. Intended use areas: Newborn, O-Vemp.

Part no.	Unit
8516225	Bag/1000 pcs (100 pcs per bag)



ABR

## Vermed Hydrogel Snap Electrode

Smaller than the rectangular version, but otherwise similar. Good soft electrode size making it very suitable for newborn screening and ABR measurements. Cloth material makes it convenient to wear. Intended use areas: Newborn, ABR.

Part no.	Unit
8102022	Bag/60 pcs (3x20 pcs)




ABR

## Vermed Tab Electrode

Rectangular tab electrode for newborn screening. 22,9x34 mm. Conductive Adhesive Hydrogel. 10x10 pcs. per pouch. Intended use areas: Newborn, O-Vemp.

Part no.	Unit
8516227	Box/1000 pcs (20x50 pcs)
8516227	bag/100 pcs



ABR

## Vermed Prewired neonatal ecg

electrode for use in pediatric and neonatal settings. 3 electrodes per pouch. square shape 22,2x22,2 mm

Part no.	Unit
8107425	Bag/300 pcs (3x100 set)



ABR

## Vermed Economy line Snap Wet Gel Electrode

Non-woven cloth electrode with conductive wet gel. Passed ISO 10993 and is Latex free. Intended use areas: Newborn, Diagnostic follow-up all ages.

Part no.	Unit
8516251	Box/1000 pcs (20x50 pcs)
8516251	Bag/50 pcs



ABR

## Vermed Allround ECG electrode

Provides very good impedance and it well suited for both newborns and adults. Can be used both as snap and tab electrode with alligator cables. Pouches come in 3 pcs, value packs available in 600 pcs. per order. Please contact for order.

Part no.	Unit
8516906	Bag/60 pcs (5x12 pcs)



8107425 Bag/300 pcs (3x100 set)

8516221 Bag/60 pcs (3x20 pcs)

8516225 Bag/1000 pcs (100 pcs per bag)

8516227 Box/1000 pcs (20x50 pcs)

8516227 bag/100 pcs

8516251 Box/1000 pcs (20x50 pcs)

8516251 Bag/50 pcs

8516906 Bag/60 pcs (5x12 pcs)





# Disposable electrodes

ABR

## Ambu® 0415M tab electrode

Medium size tab electrode with solid gel for 12-lead ECG. Presented conveniently with 10 pieces to a card. Fish tail tab design reduces electrode lift-offs. Suitable for both adults and children. Intended use areas: ABR Diagnostics of all kinds with Tab electrodes.

Part no.	Unit
8516138	Bag/100 pcs (10x10 pcs)

NEW

100 pcs



ABR

## Ambu® 2300 Blue Sensor Tab electrode

The medium size Ambu BlueSensor 2300 electrode features a skin-friendly solid gel that provides optimal signal quality recordings during 12-lead diagnostic resting ECG application. Very suitable for ABR diagnostics of all kinds. Intended use areas: ABR Diagnostics of all kinds.

Part no.	Unit
8516139	Box/1000 pcs (10x100 pcs)

NEW

1000 pcs



ABR

## Ambu® Neuroline 720

Ambu Neuroline 720 electrodes are pre-gelled and self adhesive. They are ideal for use in sleep registration where the strong, yet skin-friendly, adhesion ensures that the electrode stays in place for the entire examination, even if the patient sweats. Intended use areas: ABR Diagnostics of all kinds.

Part no.	Unit
8500060	Bag/25 pcs

25 pcs



ABR

## Ambu® Blue Sensor™ N-00-S

Foam ECG electrode for use with button cables. Wetgel.

Part no.	Unit
8013029	Bag/25 pcs

25 pcs



ABR

## Ambu® Blue Sensor™ N-10-A

Foam ECG electrode with 10 cm wire and 4 mm connector, for use with universal electrode cable. Wetgel.

Part no.	Unit
8500670	Bag/25 pcs

25 pcs



ABR

## Ambu® Blue Sensor™ P-00-S

High quality Ag/AgCl sensor – ensures excellent trace quality.

Part no.	Unit
8504783	Bag/50 pcs

50 pcs



ABR

## Ambu® Neuroline 720 – with 100 cm cable and k-connector

Ambu Neuroline 720 electrodes are pre-gelled and self adhesive. They are ideal for use in sleep registration where the strong, yet skin-friendly, adhesion ensures that the electrode stays in place for the entire examination, even if the patient sweats. Intended use areas: ABR Diagnostics of all kinds.

Part no.	Unit
8503244	Bag/10 pcs




ABR

## Ambu® Neurology Needle Electrode

This needle electrode features exceedingly low penetration-resistance levels. The needles have been coated with smooth coating which reduces friction when inserting and working in the target area.

Part no.	Unit
8519037	Box/40 pcs





# Surface electrodes tailored to your need

At Diatec we believe that it is important to choose the right surface electrode for each specific recording and patient type.

Wetgel or hydrogel electrodes?

A wetgel electrode consists of water, thickener, ionic salts, surfactants and bactericides or fungicides. Wetgel migrates easily and rapidly onto the skin surface, reducing impedance levels, which in turn will generally result in lower noise in the recording.

A hydrogel surface electrode is less adhesive when compared to wetgel, making it particularly suitable for newborns.

Chloride salt concentration

The higher the chloride salt content, the more conductive the electrode. Higher salt content (10% and above), pre-gelled surface electrodes are useful for making fast, high quality measurements of biopotentials. High conductivity electrodes generally have reduced artifact, due to the lower impedance generated between electrode and skin surface. As the chloride salt content drops, there is a corresponding reduction in electrode conductivity.

If the chloride content is below 10%, the electrode can then be employed for longer recording times, with a reduced risk of skin irritation. As a rule of thumb, wetgel has a chloride content of 10% and hydrogel 4%.

It is important to use electrodes with similar chloride salt content when testing.

The impedance of the electrode/skin junction is highly dependent on the electrolyte type and the chloride salt concentration. For example, a hydrogel electrode with 4% chloride concentration will have approximately 10x higher impedance than a wet liquid gel electrode with 10% chloride concentration, after first application to the skin.

Adhesion

An advantage of a conductive adhesive is that the active surface can be as large as the electrode itself because no additional adhesive strip is needed. The Sanibel disposable electrodes come with a range of adhesive qualities. All electrodes are designed to adhere well to skin surfaces, but the adhesion characteristics can be grouped according to the application:











- Strong adhesive
- Moderate adhesive

Strong adhesion is not recommended for newborns.

Skin preparation

For optimum electrode to skin conductivity, the skin should be lightly brushed with a gentle abrasive gel/wipe.

An alcohol wipe is not recommended to improve conductivity as this will only dry out the skin surface.

Item no.	Surface electrode	Adhesion electrode backing / Dimension (part on skin) / Material hardness / Electrode cable	Chloride conc. - Gel type	Category / Patient age	Special conditions
8107410 	<b>Sanibel™ tab electrodes</b> Packed in 20 smaller pouches of 3 electrodes (total of 60 electrodes).	Moderate / 3 x 2.2 cm / Soft tissue / Tab cables (ETSE tab surface cable kit #8011352)	Solid gel/ Hydrogel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Low chloride content, less adhesive. Especially suitable for newborns.
8107137 	<b>Sanibel™ snap electrodes</b> Pouch of 3 electrodes bundled in a box with a total of 60 electrodes.	Moderate / 3 x 2.1 cm / Soft tissue / Snap cable (ETB surface electrode kit #8011346)	Solid gel/ Hydrogel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Low chloride content, less adhesive. Especially suitable for newborns.
8107425 	<b>Sanibel™ prewired electrodes</b> 3 electrodes in each pouch with a total of 300 electrodes.	Moderate / 2 x 2 cm / Soft tissue / Prewired	Solid gel/ Hydrogel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Low chloride content, less adhesive. Especially suitable for newborns. Electrode cables are disposable, reducing cross contamination. Small soft electrode.
8108563 	<b>Sanibel™ wetgel snap electrodes</b> 12x5 pcs re-sealable pouches (total of 60 electrodes).	Strong / Ø 1 inch, 2.55 cm / Soft plastic / Snap cable (ETB surface electrode kit #8011346)	3 % - Wetgel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Strong adhesion. Low chloride content (not recommended for newborns under 1 month of age).
8109096 	<b>Sanibel™ snap electrodes</b> 60 electrodes in one resealable pouch. Repacked in 10 smaller bags of 6 pcs.	Moderate / 3.5 x 2.3 cm / Soft tissue / Snap cable (ETB surface electrode kit #8011346)	Solid gel/ Hydrogel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Low chloride content, less adhesive. Especially suitable for newborns. Large soft electrode.
8107426 	<b>Sanibel™ vinyl tab electrodes</b> 100 electrodes in one resealable pouch.	Moderate / 2.5 x 2.2 cm / Hard plastic / Tab cables (ETSE tab surface cable kit # 8011352)	Solid gel/ Hydrogel		Low cost tab electrode.
8500060 	<b>Ambu® Neuroline™ 720 side snap</b> 25 electrodes in one pouch.	Strong / 2.8 x 2 cm / Soft plastic / Snap cable (ETB surface electrode kit #8011346)	Wetgel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Strong adhesion. High chloride content (not recommended for newborns under 1 month of age).
8013029 	<b>Ambu® Blue Sensor N-00-S</b> Bag of 25 electrodes	Superior / 3 x 2 cm / Soft plastic / Button cable	Wetgel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Especially suitable for newborns.
8500670 	<b>Ambu® Blue Sensor N-10-A</b> Bag of 25 electrodes	Superior / 3 x 2 cm / foam backing	Wetgel	Specifically Developed for children	Superior adhesion to ensure optimal signal during short to medium term applications
8504783 	<b>Ambu® Blue Sensor P-00-S</b> Bag of 50 electrodes	Superior / 3.4 cm / Hard plastic / various	Wet gel	ABR, MLR, LLR, VEMP / Electrode is suitable for all ages.	Ideal for paediatric patients



ABR

### Sanibel™ ECoChG starter kit

For Eclipse™ and Bio-logic® AEP system. Includes cable, gel and 2 pcs TM electrodes.

Part no.	Unit
8013209	Set/4 pcs



ABR

### Sanibel™ ECoChG cable

For Eclipse™ and Bio-logic® AEP system. 1.5 mm/.060" female connector. For use with Sanibel™ TM electrodes.

Part no.	Unit
8500650	1 pc

### After use ...

After use, cables should be wiped with a soft cloth that has been dipped in a mild detergent solution or alcohol.

Contact between the cleaning solution and the connectors should be avoided.



ABR

### Disc electrodes

For use with Ten20™ electrode paste. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable.

Part no.	Unit
8011498	Set/3 pcs RWB
8011345	Set/5 pcs RWBYK




ABR

### Gold-plated cup electrodes

Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). Ø 10 mm gold-plated. 50 cm cable.

Part no.	Unit
8107418	Set/3 pcs RWB
8107145	Set/5 pcs RWBYK



ABR

### Alligator clip cable

For use with tab electrodes. 1.5mm/.060" female connector. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K).

Part no.	Unit
8011412	Set/3 pcs RWB
8500900	Set/3 pcs WYK
8011352	Set/5 pcs RWBYK



ABR

### Tip trode cables

For use with TTE25 gold foil tip trodes. 1.5 mm with touch proof connectors and plated clamp clips. 50 cm long.

Part no.	Unit
8011343	Set/2 pcs



ABR

### Universal electrode cable

Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable.

Non-tangling.

Part no.	Unit
8011353	Set/3 pcs RWB
8011344	Set/5 pcs RWBYK



ABR

### Standard button cable

For use with snap electrodes. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable. Non-tangling.

Part no.	Unit
8011354	Set/3 pcs RWB
8011346	Set/5 pcs RWBYK



ABR

### Alligator clip leads

For use with tab style disposable electrodes

Part no.	Description	Unit
301633	R/WT/BL, 18"	Set/3 pcs
301636	R/G/BK, 1m	Set/3 pcs
301635	R/WT/BL/G/BK, 18"	Set/5 pcs



ABR

### Pinch style cable

For use with snap electrodes. 1.5mm/ .060" female connector. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K).

Part no.	Unit
8500390	Set/3 pcs RWB
8519048	Set/3 pcs WYK
8500410	Set/5 pcs RWBYK



ABR

### Sanibel™ EP electrode jumper cable

Allows two of the preamplifier inputs to be connected, e.g. for linking mastoid electrodes for cortical testing.

Part no.	Unit
8107872	1 pc



ABR

### Reusable button snap leads 1m

Lead wires for use with snap disposable electrodes.

Part no.	Description	Unit
301601	Snap Lead, 1m R/WT/BL/G/BK	Pack/ 5pcs



ABR

### Pinch clip leads

Lead wires for use with snap disposable electrodes

Part no.	Description	Unit
301663	R/WT/BL, 18"	Set/3 pcs
301665	R/WT/BL/G/BK/, 18"	Set/5 pcs
301666	R/WT/BL/G/BK/Y, 18"	Set/6 pcs





ABR

## NuPrep™ skin prep gel, 114 g

Efficient, yet gentle skin abrasive. Chloride and acetone free.

Part no.	Unit
8500320	Box/3 tubes



ABR

## Spectra 360® electrode gel, 250 g

Salt- & chloride-free conductive electrode gel suitable for long-term applications.

Part no.	Unit
8500740	Tube/250 g



ABR

## Lectron® II electrode gel, 250 ml

Highly conductive, chloride-free electrode gel. Recommended for use with Sanibel™ TM electrodes for ECoG.

Part no.	Unit
8102626	Tube/250 ml



ABR

## Everi™ abrasive paste, 160 g

Reduces the skin impedance during EEG, EP, ECG and Holter Stress test examination.

Part no.	Unit
8012952	Bottle/160 g



ABR

## SignaGel®, 250 g

Hypoallergenic, bacteriostatic conductive saline electrode gel.

Part no.	Unit
8102628	Tube/250 g



ABR

## Sonavelle® electrode gel, 250 ml

Formaldehyde-free and glycerol-free.

Part no.	Unit
8500380	Bottle/250 ml



ABR

## Ten20™ conductive EEG paste, 114 g/228 g

Adhesive gel for disc electrodes. Available in pack of 3x4 oz/114 g tubes, 3x8 oz/228 g jar.

Part no.	Unit
8012924	Box/3 tubes
8500330	Box/3 jars



ABR

## 3M™ Red Dot™ trace prep, 5 m

One-step skin preparation. 18 mm x 5 m / 0.70" x 196".

Part no.	Unit
8500030	Dispenser/5 m



## Patient preparation



12 rolls



ABR

### 3M™ Transpore™ Surgical tape, 9.14 m

Transparent, latex-free, water-resistant and porous for breathability.  
25 mm x 9.14 m / 1" x 10 LY.

Part no.	Unit
8500750	Box/12 rolls



12 rolls



ABR

### 3M™ Micropore™ Surgical tape, 9.14 m

Latex-free and hypoallergenic. Holds well on damp skin.  
25 mm x 9.14 m / 1" x 10 LY.

Part no.	Unit
8500760	Box/12 rolls



ABR

OAE

### Gauze swabs

Non-sterile, 7.5 x 7.5 cm and 8 layered.

Part no.	Unit
8500040	Pack/100 pcs.



10 packs



ABR

OAE

### Q-tips

Cotton tipped applicators for patient preparation. Non-sterile.  
Each pack contains 100 pcs.

Part no.	Unit
8500810	Box/10 packs



ABR

OAE

### Small Alcohol Pads

70% isopropyl alcohol. Breimer swabs  
30 x 110mm

Part no.	Unit
9581200	Box/200 pcs.



500 pcs



### 3" Disposable ear cushion cover

Replaces disinfectants. Headphone Covers for classic audiometric headphones.

Part no.	Unit
8107419	Box/500 pcs



### Probe cleaning tool

For Interacoustics Titan™/AT235™/AA222™/Sera™/MT10™/Amplivox Otowaves/MAICO Mi XX + touchTym™ and GSI TymStar™ Pro.

Part no.	Unit
8011351	1 pc



### Disinfectant detergent, 1000 ml

Alcohol-free, ready to use solution for a quick and extensive disinfection of non-invasive medical devices.

Part no.	Unit
8501706	Bottle, 1000 ml

100 pcs



### 4.5" Disposable ear cushion cover

Replaces disinfectants. Headphone Covers for noise enclosure audiometric headphones. For HDA200/DD450/DD65v2/Peltor headsets.

Part no.	Unit
8513404	Bag/100 pcs



### Keyboard protective foil

Place a protective foil over your keyboard to minimize contamination with grime, bacteria and germs, or liquid.

Part no.	Unit
8500460	1 pc



### Clear safety glasses

Part no.	Unit
9000564	1 pc



ABR IMP OAE

## Bridge & Implant floss

Used to clean probe tips.

Part no.	Unit
8003566	Box/30 pcs



ABR IMP OAE

## 3-in-1 floss

Used to clean probe tips.

Part no.	Unit
8003565	Box/100 pcs



ABR IMP OAE

## Spiraflex Cleaners

Used to clean probe tips.

Part no.	Unit
8902124100	Bag/30 pcs




## Large alcohol swabs

70% isopropyl alcohol, 2% CHG. PDI swabs.  
150 x 162 mm

Part no.	Unit
900331	Box/200 pcs

How to clean  
probe tips!







## PDI Alcohol wipes canister

70% isopropyl alcohol

Part no.	Unit
900325	200 sheets



## Canister refill

70% isopropyl alcohol. Suitable for use with 900325

Part no.	Unit
900328	200 sheets



## Spirigel sanitising hand gel, pump

Kills germs without water, air dries in seconds, added moisturiser.

Part no.	Unit
9045940	500ml



## PDI nail polish remover pads

Part no.	Unit
9005565	Box/100 pads



## PDI Detergent wipes, bucket

Wipe size: 300 x 245mm

Part no.	Unit
9000175	Bucket/225 pcs



## PDI Detergent wipes, refill

Solo lid bag. Wipe size: 300 x 245mm.

Part no.	Unit
8000210	Bag/225 pcs

# Infection control

## A clean clinic is also part of your professional image

Infection control is so much more than cleaning probe tips, protecting your personal computer's keyboard, and maintaining high standards within hygiene using e.g. disinfectant detergent for a quick and extensive disinfection of non-invasive medical devices.

When it comes to infection control, preventing germs, dirt and dust from spreading via other sources, is also of critical importance – both to the safety of your patients and yourself! Gloves, shoe covers, facemasks, and aprons are just some of the “accessories”, which can help you succeed in this.

Maintain high, professional standards and make your clinic an ever-safe environment for your patients, making the wearing of gloves, facemasks and other body protection wear a natural part of the “patient interaction” process.



## Bouffant Caps

Perform high Hygiene practices in front of your patients or clients.

Part no.	Unit
8519350	Bag/250 pcs (21 Inches)



## Disposable facemasks with loops

Perform high Hygiene practices in front of your patients or clients.

Part no.	Unit
8515555	Pack/50 pcs

## Disposable apron

Keep a professional appearance also after servicing a hearing instrument.

Part no.	Unit
8519348	Box/100 pcs



## Promed nitrile soft examination gloves

Latex and powder-free and provide excellent protection. Our Nitrile gloves are available in many sizes and will accommodate most hands..

Part no.	Description	Unit
8515551	Small	Box/100 pcs
8515552	Medium	Box/100 pcs
8515553	Large	Box/100 pcs
8515554	Extra Large	Box/100 pcs \$20

## Disposable Shoe Covers

Perform high Hygiene practices in front of your patients or clients.

Part no.	Unit
8519349	Box/100 pcs (Regular)





24 pcs

BAL

### Disposable VNG goggle foam pads

For Interacoustics® VNG mask.

Part no.	Unit
8011527	Box/24 pcs



New

### Disposable specula for Difra Airstar

Part no.	Description	Unit
8011457	2.5 mm spec. (green)	bag/100 pcs
8011456	4 mm spec (blue)	bag/100 pcs



New

### Paper Kidney tray

Part no.	Description
8520139	10 pcs



34 pcs


BAL

### Disposable ear specula

For AirFx™.



Part no.	Description	Unit
8013013	2.75 mm	Bag/34 pcs
8013014	4.25 mm	Bag/34 pcs



50 pcs

BAL

### Water irrigator tips

For use with Interacoustics® Aqua Stim™.



Part no.	Unit
8013018	Bag/50 pcs



New

AUD IMP

### Accessories for Otoscopy

Stylish wall dispenser for ear funnels. Can contain 6 different sizes.

Part no.	Description
8515678	Specula wall dispenser



No sharp edges

100 pcs

AUD BAL IMP

### Video Otoscopes specula

Fits the new Interacoustics® Viot and OriVision otoscope.

Part no.	Description	Unit
8515083	KaWe 2.5 mm speculum	Bag/100 pcs
8511728	KaWe 4 mm speculum	Bag/100 pcs



### Disposable ear specula

For Welch Allyn.

Part no.	Description	Unit
8107750	2.5mm	Bag/250 pcs
8108252	4.0mm	Bag/250 pcs

Box of 1000 pcs available



### Disposable ear specula

For Heine Otoscope. Suitable for adult and children.

Part no.	Description	Unit
8450102907	2.5 mm	Bag/50 pcs
8450102806	4.0mm	Bag/50 pcs



No sharp edges

50 pcs

AUD BAL IMP

### Disposable speculas for video otoscope

Fits the AURICAL™ OTOcam 300, Interacoustics® VIOT v1 and XION™ video otoscopes.

Part no.	Unit
8500470	Bag/50 pcs



### Bulb for Heine Mini Otoscope

2.5v Halogen replacement bulb 110

Part no.	Unit
8900204705	1 pc



# Heine Mini Otoscope




BAL

## Heine Mini Otoscope

Heine Mini 3000 set. Includes 5 speculas.

Part no.	Unit
8900400105	Box Set

New



AUD

IMP

## Luxamed Otoscope Auris LED 2.5V

Otoscope from German Luxamed. LED light to ensure longevity and great light source, dimable. Fits Heine/KaWe compatible speculums.

Part no.	Unit
8519433	1 pc

## Luxamed Otoscope Auris LED 3.7V USB/PC Charged

Otoscope from German Luxamed. LED light to ensure longevity and great light source. USB charging. Fits Heine/KaWe compatible speculums.

Part no.	Unit
8519434	1 pc



**Round nose scissors**

Part no.	Unit
9100234	1 pc



**Mini pliers**

Part no.	Unit
9100235	1 pc



**Tube expander**

Part no.	Unit
930421	1 pc



**Mini wire cutters**

Part no.	Unit
9100236	1 pc



**Cleaning tool**

Part no.	Unit
8902126000	Bag/20 pcs



**Tweezers**

Part no.	Unit
9100237	1 pc



OAE

## Thermal paper for ERO•SCAN™/OtoRead™, 25 m

Width: 57 mm

Inner diameter: 12 mm

Outer diameter: 45 mm

Part no.	Unit
8500200	1 roll/25 m



ABR IMP OAE

## API300 thermal paper for Titan™, 10 m

Width: 57 mm

Inner diameter: 12 mm

Outer diameter: 35 mm

Part no.	Unit
8500230	1 roll/10 m



ABR IMP OAE

## Sanibel™ MPT-III thermal paper, 20 m

For Sanibel™ MPT-III printer which is supplied with Interacoustics® AT235™.

Width: 79 mm

Inner diameter: 12 mm

Outer diameter: 45 mm

Part no.	Unit
8500610	1 roll/20 m



IMP OAE

## Thermal paper for Sanibel™ MPT-II printer, 20 m

For the Sanibel™ MPT-II printer which is supplied with MAICO® ERO•SCAN™, Grason-Stadler® Corti™ and Interacoustics® OtoRead™.

Width: 57 mm

Inner diameter: 9 mm

Outer diameter: 40 mm

Part no.	Unit
8029305	1 roll/20 m



IMP

## TPR10™ thermal paper for MTP10™, 25 m

Width: 79 mm

Inner diameter: 12 mm

Outer diameter: 45 mm

Part no.	Unit
8500220	1 roll/25 m



AUD IMP

## TPR26™ thermal paper for Interacoustics®, 25 m


Fits most popular audiometers and tympanometers.

Width: 112 mm

Inner diameter: 12 mm

Outer diameter: 45 mm

Part no.	Unit
8500210	1 roll/25 m



IMP

## Adhesive thermal paper LPQ80, 195 labels

For GSI TympStar™ Pro.

Width: 82 mm

Inner diameter: 25 mm

Outer diameter: 88 mm

Part no.	Unit
8509894	1 roll/195 labels



IMP

## Thermal paper LPQ80, 80 m


For GSI TympStar™ Pro.

Width: 80 mm

Inner diameter: 13 mm

Outer diameter: 80 mm

Part no.	Unit
8509893	1 roll/80 m



IMP

## Thermal paper for GSI70

For GSI 70.

Width: 109.5 mm

Part no.	Unit
8500240	1 roll



IMP

## f chart AF12

100 pcs

Part no.	Unit
8012858	1 pad



For Scriber Printers.

Width: 56 mm

Length: 60 mm

Part no.	Unit
8107448	1 roll/ 40 labels



## For OtoRead™ PRT MPT-II printer, 10 m

Part no.	Unit
8013009	1 roll/10 m





## Paper



IMP

### Thermal paper for GSI

117 Thick, Adhesive,  
9m long. For GSI Tymptstar

Part no.	Unit
8012937	1 roll



IMP

### GSI Thermal paper

For GSI39  
Width: 111 x 51 x 12 mm

Part no.	Unit
8012939	1 roll



IMP

### GSI Thermal paper with label

For GSI37  
Width: 58 mm

Part no.	Unit
8012950	1 roll



IMP

### Audioscan printer paper

For VFI, Axiom and RM500SL  
Width: 80 mm

Part no.	Unit
RE201-25	1 roll

## Wax removal accessories



### Bionix Safe Ear Currettes

Safe Ear Currettes allows safe and gently removal of different types of cerumen. Single use ear currettes. Ensure your patients' safety by using these affordable ear currettes. Available as variety pack or single size items. Safe Ear Currettes allow healthcare providers to safely and gently remove cerumen with styles to match different types of wax and personal preference.

Part no.	Description	Unit
8515585	Green MicroLoop®	50 pcs/box
8515586	Orange ControLoop®	50 pcs/box
8515587	Red AngleLoop®	50 pcs/box
8519359	Purple VersaLoop®	50 pcs/box
8519360	Blue InfantScoop®	50 pcs/box
8519372	Yellow CeraSpoon®	50 pcs/box
8519361	White FlexLoop®	50 pcs/box
8519362	Variety Kit, (10 each style & 15 white)	



### Lighted Foreceps

Bionix Lighted Foreceps combine illumination and magnification to help healthcare providers take less time to see and safely remove ear canal and nasal passage obstructions.

Part no.	Description
8519392	10 foreceps, 1 light and 1 magnifying lens



### Bionix Lighted Ear Currettes

Box of 50 pcs. Use the Bionix lighted ear currettes for safe removal of earwax. The Bionix currettes comes in a variety of different shapes for you to option for your optimal choice in a situation. Bionix Lighted Ear Currettes combine illumination and magnification to help healthcare providers take less time to see and safely remove cerumen from the ear canal.

All versions including 1 Light Source and 1 magnifying lens.

Part no.	Description
8519365	Lighted AngleLoop®, 3mm diameter
8519364	Lighted MicroLoop®, 4mm diameter
8519366	Lighted WaveCurette™, 2.5mm wide tip
8519367	Lighted FlexLoop®, 4mm wide tip
8519368	Lighted CeraSpoon®, 4mm wide tip
8519369	Lighted InfantScoop®, 2mm wide tip
8519373	Lighted VersaLoop®, 3mm wide tip
8519370	Lighted CeraPik™, 2.5mm concave tip. Great for hard, impacted cerumen and foreign bodies.



### Lighted Ear Curette Variety Pack

Part no.	Description
8519363	with 1 Light source and magnifying lens includes all tip styles except CeraPik™

## Wax removal accessories



### Earlight for wax removal

Part no.	Description	Unit
8515592	LED Earlight Straight Tip	Bag/4 pcs
8515594	LED Earlight Pros Tip	Bag/4 pcs
8515595	LED Earlight Scoop Tip	Bag/4 pcs
8515596	LED Cerumen Management TIP Stainless Steel Loop	Bag/2 pcs
8515597	LED Earlight (with 2 batteries)	



### OtoClear Irrigation Tips

Bionix OtoClear Ear Irrigation helps provide safe and comfortable ear irrigation without the backslash, at half the cost of a self-constructed device.

Part no.	Unit
8519375	Box 40 pcs

### Clean ears spray, 30ml

Dissolves earwax and irrigates ear canal. Simple easy-to-use spray which is suitable for all ages.

Part no.	Unit
9001258	1 pc



### LED Earlight

Part no.	Unit
9125718	1 set



### OtoClear Spray wash Kit

OtoClear Ear Irrigation system directs water to the ear canal wall gently with three divergent streams, instead of directly at the Tympanic Membrane

Part no.	Description
8519371	box 20 irrigator tips, 1 470ml spray wash bottle, 1x ear basin



### OtoClear® ear basin

Part no.	Unit
8519390	Box/3 pcs



### Earlight with Scaled Top

Part no.	Description	
8519357	Luxamed Ear light black, acrylic scaled top	\$5.00
8519358	Luxamed Acrylic Top with Scale	\$45.00

New

## Wax removal accessories



### Micro Suction Products

Otopront VAC18 is a transportable suction unit, specific designed for ear wax removal or other temporary suction purposes within hospitals and clinics.

Part no.	Description
8515194	Otopront VAC 18 suction unit



### Accessories for Otopront Suction units

To optimize the use of the Otopront VAC18. Please use the appliances below per demand.

Part no.	Description
8517966	Suction container system (to be added the VAC 18)
8515200	Refill bags for disposables suction system 60 pcs / box, Value box of refill bags
8515372	Refill bags for disposables suction system, 10 pcs
8520345	Silicon suction hose medical grade 50m
8515195	Air filter VAC 18
8515197	Internal filter VAC 18 / VAC 28
8515199	Lid for suction container VAC 18 / VAC 28


New



### EarWay Pro - Wax Removal Helix

Very easy and safe way to remove ear wax. The helix grabs around the ear wax lump and pulls it out in a comfortable way.

Part no.	Description	Unit
8518025	Pro M, 5.4 mm, Purple,	25 units
8518024	Pro S, 4.2 mm, Blue,	25 units







AUD

## USB cables A to B

Part no.	Description	Unit
8011241	2 m/6.56 feet	1 pc
8011243	5 m/16.40 feet	1 pc



AUD

## HDA200 replacement cable, symmetric

For HDA200 and DD65 headphones. 2 m cable and 6.3 mm/ 30° jack plug.

Part no.	Unit
8004172	1 pc



AUD

## TDH39™/DD45™ replacement cable, symmetric

For TDH39™ and DD45™. 30° jack plug.

Part no.	Description	Unit
8010822	2 m	1 pc
8507204	5 m	1 pc



AUD

IMP

## TDH39™/DD45™ contra headset replacement cable

For TDH39™ and DD45™. Single mono. 30° jack plug.

Part no.	Unit
8010826	1 pc



AUD

## B71 Bone conductor replacement cable

2 m cable and 6.3 mm/ 30° jack plug.

Part no.	Unit
8011094	1 pc



AUD

## IP30 Contra Insert earphone cable

Cable length: 2 m. Plug type: 30°, 1x6.3 mm gold-plated mono jack.

Part no.	Unit
8107358	1 pc



## Parts for audiometry



AUD

### RadioEar Foldable audiometric headband

Black, foldable.

Part no.	Unit
8010882	1 set



AUD

### RadioEar Headband with 2 yokes

Black.

Part no.	Unit
8010840	1 set



AUD

### Peltor® headset/earmuff replacement hygiene kit

For DD450, HDA200 and DD65.

Part no.	Unit
8500270	1 set



AUD

### RadioEar DD45 Audiometric headset

Headband: HB7. Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8010954	1 set



AUD

### RadioEar DD45 Adult audiometric headset

Headband: HB3045. Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8105268	1 set

### RadioEar DD45 Adult audiometric headset

Headband: HBA, foldable. Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8106340	1 set



AUD

### RadioEar DD65 v2 Passive attenuation audiometric headset

Cable length: 2 m. Plug type: Straight, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8510117	1 set



AUD

### RadioEar DD45 transducer

Audiometric earphones engineered to match the performance and function of the most popular earphones available for audiological testing. 10 ohm receiver. Easy installation. Reliable. RoHS compatible.

Part no.	Unit
8507268	1 set

AUD

### RadioEar DD45 Contra headband

Tension: 4,5 N ±0.5N

Part no.	Unit
8507211	1 set



Medical grade  
silicone rubber  
Biocompatible



AUD

### Sanibel™ TDH39/DD45 earphone cushion red & blue, model 51

For TDH39, 49, 50 & DD45 earphones.

Part no.	Unit
8107611	1 Set



Medical grade  
silicone rubber  
Biocompatible



AUD

### DD 51 cushion

For DD45™ and TDH39™ headsets.

Part no.	Unit
8010894	1 pc







AUD

## Foam pad for contra headset

Compatible with e.g. TDH39c™ and DD45c™.

Part no.	Unit
8507292	Bag/10 pcs



AUD HIF ABR

## RadioEar replacement pads (3288-1)

For RadioEar Spring Steel headband.

Part no.	Unit
8504777	Bag/10 pcs



AUD

## Ear cushion HZP 09 pair

Ear cushion foam (1 pair) suitable for: PC 20, PC 25, PC 30, PC 31, PC 35, PC 36, PC 130, PC 131, PC 135, PC 136, PC 140, PC 141, PC 145, PC 146.

Part no.	Unit
8104416	1 pair



AUD

## Headband cushion

For RadioEar HBA foldable audiometric headband.

Part no.	Unit
8010886	1 pc



AUD

## RadioEar B71W Bone transducer headset

Ohm: 10 ohm. Cable length: 79"/200 cm molded cord. Plug type: 30°, 1x6.3 mm gold-plated mono jack. Headband: P3333

Part no.	Unit
8510370	1 pc



AUD HIF ABR

## RadioEar Spring Steel headband (P3333)

Part no.	Unit
8507213	1 pc



AUD

## Headband cushion

For Sennheiser HDA200, RadioEar DD450, DD65v2 and Peltor audiometric headbands.

Part no.	Unit
8000770	1 pc



10 pcs

AUD HIF ABR

## RadioEar AMBAND bone conductor headband (P1188)

For adult, children and infants.

Part no.	Unit
8507214	Bag/10 pcs




New

## Replacement cable for bone conductor B71/B81 (with safety fixation point)

Part no.	Unit
8107612	1 pc





New

### Apex 2 Inlet Filters

Part no.	Unit
209127A	Bag/10 pcs



New

### OAE Probe body nozzle

Part no.	Unit
580-PRBNZL	1 pc



New

### Ultra-Vac Needle pack

Part no.	Unit
8111306	Bag/6 pcs



New

### Casella Apex Inlet Foam Filter

Part no.	Unit
P104122	Bag/10 pcs



New

### Ultra-Vac Filter

Part no.	Unit
8111312	1 pc



HIF

### Ear coupler sealing wax

Use to create a seal between coupler and hearing instrument.

Part no.	Unit
8500080	1 pc



ABR AUD OAE

### Alkaline battery 1.5V AA

Part no.	Unit
8011226	1 pc



### Clothing clip kit

Part no.	Unit
8103992	1 pc



### Eartip storage cabinet

12 Small drawers, 3 Medium drawers, 1 Large drawer. Can be hung on the wall. H 28.2 cm x W 30.6 cm x L 15 cm

Part no.	Unit
8100665	1 pc



### Battery tester

Part no.	Unit
8901014804	1 pc



IMP HIF

### Lithium battery 3.7V

For Interacoustics® Callisto™ and Titan™.

Part no.	Unit
8009978	1 pc



5 pcs

AUD IMP

### EARNimals®

Cartoon faces that are placed on the otoscope. Make ear examinations less scary for kids.

Part no.	Unit
8011501	Bag/5 pcs (assorted)




### CR2032 Battery

Lithium Cell

Part no.	Unit
9000132	1 pc



To order, please send us an e-mail  
[info@diatec-diagnostics.co.nz](mailto:info@diatec-diagnostics.co.nz) or  
contact us on 09 279 8833



Diatec New Zealand • [info@diatec-diagnostics.co.nz](mailto:info@diatec-diagnostics.co.nz)  
Diatec is a company within the William Demant Group

**Disclaimer**

This catalogue has been published for worldwide distribution. Some parts described may NOT BE AVAILABLE in your area. Your local distributor will be pleased to give you details. For advice on compatibility and instructions of use, please consult your operation manual or local distributor. Sanibel Supply A/S assumes no legal liability or responsibility for the accuracy, typing errors or mistakes, correctness or completeness of any information in this flyer. The information in this flyer was correct to the best of our knowledge at the time of printing. Sanibel Supply A/S reserves the right to make changes at any time in specification of products described in this flyer without notice or obligation. MADSEN AccuScreen® and MADSEN OTOflex® are registered trademarks of Otometrics A/S. Echo-Screen®, AuDx®, Navigator® Pro, and ABAer® are registered trademarks of Natus Medical Incorporated. Sanibel™ is a registered trademark of Sanibel Supply®. All other Instrument brands or model names are copyrighted or trademarked by their respective owners.

