
Diatec Consumables Catalogue 2021





Contents



Page 4 - 7

Disposable ear tips



Probe tips



Page 10 - 12

Supplies for insert earphones



Page 13 - 14

ABR & newborn



Page 15

Rem probe tubes for HA fitting



Page 16 - 18

Hearing Aid Fitting & supplies



Page 19 - 20

Hearing Aid Maintenance & hearing care



Page 21 - 24

Impression materials & supplies



Page 25 - 29

Disposable electrodes



Page 30 - 31

Electrode cables



Page 32 - 33

Patient preparation



Page 34 - 36

Infection control



Page 37

Personal Protective Equipment



Balance supplies



Page 39 - 40

Otoscopy Supplies



Page 41

Clinical tools



Page 42 - 43

Paper



Page 44 - 50

Wax removal accessories



Page 51 - 55

Parts for audiometry



Page 56

Miscellaneous

List of abbreviations



Hearing Instrument Fitting

Impedance

Auditory Brainstem Response



Audiometry

Otoacoustic Emissions



Single patient use only!



Product items are not made with natural rubber latex

Selecting **the correct** ear tip





Selecting the correct size of ear tip takes practice

In general, flanged or mushroom-shaped ear tips will be used for most testing (Impedance, DPOAE, TEOAE and Automated ABR).

Sometimes, an umbrella-shaped ear tip will be used for performing a quick impedance screening.

The table can be used as a guide to help you select the appropriate ear tip size.

Disposable ear tips

Disposable Sanibel™ ADI ear tips

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip: Amplivox: all OtoWave versions. Interacoustics: Titan, AT235, AA220, Eclipse, MT10, Sera. Maico: EasyTymp, MI 22, MI 23, MI 24, MI 26, MI 34, MI 44, Easyscreen



ADI Assortment Box

Part no.	Description U	nit
8511445	Various sizes for Tympano	metry and OAE
8501572	Screening ear tip kit, vario screening	us sizes for newborn



Flanged infant ear tips

Part no.	Description	Unit
8502838	Pediatric (green)	Bag/100 pcs
8012964	3-5 mm (red)	Bag/100 pcs
8012966	4-7 mm (blue)	Bag/100 pcs
8012968	5-8 mm (yellow)	Bag/100 pcs









Umbrella ear tips

Part no.	Description	Unit
8012986	7.5 mm (green)	Bag/100 pcs
8012988	9 mm (red)	Bag/100 pcs
8012990	10 mm (blue)	Bag/100 pcs
8012992	11 mm (yellow)	Bag/100 pcs
8012994	13 mm (green)	Bag/100 pcs
8012996	15 mm (red)	Bag/100 pcs
8012998	19 mm (blue)	Bag/100 pcs
8013000	22 mm (yellow)	Bag/100 pcs











IA series ear tips

For use with many acoustic immittance instruments and otoacoustic emission systems, including: GSI: TympStar Pro, Allegro, Novus and Corti. Interacoustics (clear probe): AT235, AA222 and Titan and Otoread II, MT10. Maico: Touchtymp Racecar, MI24, MI26, MI34, MI36, Easytymp. Madsen/Otometrics: Zodiac (new), Accuscreen. Path Medical: all Sentiero products,





Infant ear tips

Part no.	Description	Unit
8013006	3 mm (yellow)	Bag/100 pcs
8002035	4 mm flanged (red)	Baa/100 pcs





Mushroom ear tips

Part no.	Description	Unit
8013002	7 mm (blue)	Bag/100 pcs
8013004	8 mm (yellow)	Bag/100 pcs
8012970	9 mm (green)	Bag/100 pcs
8012972	10 mm (red)	Bag/100 pcs
8012974	11 mm (blue)	Bag/100 pcs
8012976	12 mm (yellow)	Bag/100 pcs
8012978	13 mm (green)	Bag/100 pcs
8012980	14 mm (red)	Bag/100 pcs
8012982	15 mm (blue)	Bag/100 pcs
8012984	19 mm (yellow)	Bag/100 pcs





Mushroom ear tips

Part no.	Description	Unit
8520348	3-6 mm (green)	Bag/100 pcs
8520349	6 mm (yelllow)	Bag/100 pcs
8520350	7 mm (red)	Bag/100 pcs
8520351	8 mm (blue)	Bag/100 pcs
8520352	9 mm (green)	Bag/100 pcs
8520353	10 mm (yellow)	Bag/100 pcs
8520354	11 mm (red)	Bag/100 pcs
8520355	12 mm (blue)	Bag/100 pcs
8520356	13 mm (green)	Bag/100 pcs
8520357	14 mm (yellow)	Bag/100 pcs
8520358	16 mm (red)	Bag/100 pcs
8520359	18.5 mm (blue)	Bag/100 pcs





IT Select Series Ear Tips

For use with Inventis: Flute, Clarinet, Viola. Resonance: r16m, r26m, r36m, r15c, r25c, r35c











IT Select

Part no.	Description	Unit
8522402	8mm (green)	Bag/100 pcs
8522403	10mm (red)	Bag/100 pcs
8522404	12mm (blue)	Bag/100 pcs
8522405	14mm (green)	Bag/100 pcs
8522406	16mm (red)	Bag/100 pcs







AKR Mushroom ear tips

Use with the following acoustic immittance instruments and otoacoustic emission systems with BLACK probe tip from Interacoustics and Maico: AT235, AA220, AA220, Eclipse, MT10, Otoread (I), RaceCar MI 22, MI 23, MI 24, MI 26, MI 34, MI 44. Madsen/Otometrics: Capella 2, Zodiac 901, Biologic AuDX Pro (+Flex). Neurosoft: Neuro-Audio, Oscilla: TSM 300/400 + T800 series. All GSI tympanometers from GSI 20 - 39 + Tympstar.



Mushroom ear tips

Part no.	Description	Unit
8518409	AKR Assortment Box, various sizes	<u> </u>
8500570	3-6 mm (yellow)	bag/100 pcs
8107510	4 mm (red)	bag/100 pcs
8107511	5 mm (blue)	bag/100 pcs
8500560	6 mm (green)	bag/100 pcs
8500550	7 mm (yellow)	bag/100 pcs
8500540	8 mm (red)	bag/100 pcs
8500530	9 mm (blue)	bag/100 pcs
8500520	10 mm (green)	bag/100 pcs
8500510	11 mm (yellow)	bag/100 pcs
8500500	12 mm (red)	bag/100 pcs
8500490	13 mm (blue)	bag/100 pcs
8107520	14 mm (green)	bag/100 pcs
8500480	16 mm (yellow)	bag/100 pcs









Snap probe accessories

Accessories specifically designed to work with SnapPROBE*. For use with **Interacoustics** Sera, **GrasonStadler** Novus, **Maico** EasyScreen.

Part no.	Description	Unit
8522758	Ear tip Infant A (Conical)	Bag/100 pcs
8522759	Ear tip Infant A (Conical)	Bag/25 pcs
8522760	Ear tip Infant C (Narrow)	Bag/100 pcs
8522761	Ear tip Infant C (Narrow)	Bag/25 pcs
8516436	Probe Waxguard	1 pc
8517889	Probe Waxguard	Set 10 pcs
8522461	Assortment of accessories including	
	Infant tip A and C and wax	kguard



Sanibel™ AZE mushroom ear tips

For use with old MT10 and the following systems with black probe tip: AT235, AA220, MI22, MI23, MI24, MI26, MI34, MI44, and old SPL60 (4mm)

Part no.	Description	Unit
8107622	6 mm (violet)	Bag/100 pcs
8107623	7 mm (yellow)	Bag/100 pcs
8107625	10 mm (blue)	Bag/100 pcs
8107626	11 mm (white)	Bag/100 pcs
8107629	18 mm (blue)	Bag/100 pcs









AOT ear tips

For use with Madsen/Otometrics AccuScreen*, Biologic Bio-logic* AuDX*. Pathmedical: all Sentiero products

Part no.	Description	Unit
8106367	3.7 mm (yellow)	Bag/50 pcs
8106368	4.0 mm (light green)	Bag/50 pcs
<u>8</u> 106369	4.5 mm (light pink)	Bag/50 pcs
8106370	5.0 mm (light blue)	Bag/50 pcs
8106371	4-7 mm (transparent)	Bag/50 pcs





DISPOSABLE SANIBEL™ AOT and ANA ear tips



ANA ear tips for echo-Screen®

For use with Natus: echo-Screen*. Madsen/ Otometrics: AccuScreen* Classic. Interacoustics: Eclipse with pink probe tip

Part no.	Description	Unit
8106363	4.0 mm (blue)	Bag/50 pcs
8106362	4.5 mm (red)	Bag/50 pcs
8106364	5.0 mm (green)	Bag/50 pcs
8106366	Infant tree (yellow)	Bag/50 pcs
8106365	Mushroom (green)	Bag/50 pcs







Grason & Associates ot & MO select series ear tips



OT select ear tips w/ triangle holes

Ear tips with triangle holes for use with Otodynamics UGD, SND, SGD probes

Part no.	Description	Unit
8102644	4.8 mm (blue)	Bag/100 pcs
8102645	6 mm (green)	Bag/100 pcs
8102646	8 mm (red)	Bag/100 pcs
8102647	10 mm (blue)	Bag/100 pcs
8102648	12 mm (green)	Bag/100 pcs
8102649	14 mm (red)	Bag/100 pcs
8102650	16 mm (blue)	Bag/100 pcs







OT select ear tips w/two probe holes

Ear tips with two probe holes for use with Otodynamics ECP, UGS, SNS, SGS probes.

Part no.	Description	Unit
8102632	3.5 mm (blue)	Bag/100 pcs
8102633	4.5 mm (blue)	Bag/100 pcs
8102634	5.5 mm (blue)	Bag/100 pcs
8102635	6.5 mm (blue)	Bag/100 pcs
8102636	4.8 mm (blue)	Bag/100 pcs
8102637	6 mm (green)	Bag/100 pcs
8102638	8 mm (red)	Bag/100 pcs
8102639	10 mm (blue)	Bag/100 pcs
8102640	12 mm (green)	Bag/100 pcs
8102641	14 mm (red)	Bag/100 pcs
8102642	16 mm (blue)	Bag/100 pcs







MO select ear tips for Madsen OTOflex

For MADSEN® OTOflex™.

Part no.	Description	Unit
8102653	3-6 mm (yellow)	Bag/100 pcs
8102654	6 mm (red)	Bag/100 pcs
8102655	7 mm (blue)	Bag/100 pcs
8102656	8 mm (green)	Bag/100 pcs
8102657	9 mm (yellow)	Bag/100 pcs
8102658	10 mm (red)	Bag/100 pcs
8102659	11 mm (blue)	Bag/100 pcs
8102660	12 mm (green)	Bag/100 pcs
8102661	13 mm (yellow)	Bag/100 pcs
8102662	14 mm (red)	Bag/100 pcs
8102663	16 mm (blue)	Bag/100 pcs
8102664	18.5 mm (green)	Bag/100 pcs





Otodynamics Ear Tips

For UGD or SND probes

Part no.	Description	Unit
T-R4.8S	4.8mm (blue)	Bag/100 pcs
T-R5.8B	5.8mm (blue)	Bag/100 pcs
T-R6.8B	6.8 mm (blue)	Bag/100 pcs



Otodynamics Ear Tips

Eartips for measuring DPOAE and TEOAE. For use with Otodynamics ILO-OAE instruments with DPOAE and TEOAE probes.

Part no.	Description	Unit
T-R7M	7mm (blue)	Bag/100 pcs
T-R8M	8mm (blue)	Bag/100 pcs
T-R9M	9mm (blue)	Bag/100 pcs
T-R11M	11mm (blue)	Bag/100 pcs
T-R13M	13mm (blue)	Bag/100 pcs

Ear tips for OAE screeners

For use with Biologic AuDx, Welch Allyn OAE screener and Natus echo-Screen* III.



OAE ear tips, AuDx

Part no.	Description	Unit
8107885	AuDX Preemie 3 mm	Bag/100 pcs
8107889	AuDX Tree 3.5-5mm	Bag/100 pcs
8107887	AuDX Infant 6 mm	Bag/100 pcs
8102681	Adult OAE foam	Bag/100 pcs
8102682	Pediatric OAE foam	Bag/100 pcs
8103123	Jumbo OAE foam	Bag/100 pcs





Grason & Associates mixed ear tips



TS211

Smaller offices needing the full size selection of 4-16mm may prefer the TS211 with 10 of each KR ear tip, sizes 4-16mm.



TS291

Smaller offices needing the full size selection of 3-14mm may prefer the TS291 with 10 of each GSI mushroom ear tip, sizes 3-14mm.



DP Probe ear tips

Large foam tips for adult ears.

Part no.	Unit
T-R100L	Bag/100 pcs



Using a new ear tip is optimal for obtaining a seal when inserting it into the patients' ear. A multi-used ear tip may cause inaccuracies during testing.



Clinical Audiology Kit

Suitable for students, graduates and new clinics

Part no.	Description	Unit	
1100141		1 kit	
	Includes:		
	1x Heine Mini 3000	Otoscope	
	1x LED Earlight		
	50x Heine Disposable Specula 4.0mm		
	50x Heine Disposable Specula 2.5mm		
	1x Stetoset		
	1x Combo 3 in 1 tool		
	1x Vent cleaner		
	5x Sanibel Universal Type A Probe tube		
	20x Alcohol Pads		
	20x 10mm Green impedance tips		
	20x 11mm Yellow impedance tips		
	20x 12mm Red impedance tips		
	20x 13mm Blue imepdance tips		
	20x 14mm Green ir	mpedance tips	
	20x 3A foam tips		
	20x 3B foam tips		
	20x Disposable ea	r cushion covers	

2x Clear plastic containers



Eartip storage cabinet

12 Small drawers, 3 Medium drawers, 1 Large drawer. Can be hung on the wall. H 28.2 cm x W 30.6 cm x L 15 cm

Part no.	Unit	
8100665	1 pc	







Probe tip for Sanibel™ ADI™ ear tips

Fits Otowave and new MT10™.

Part no.	Unit	
8002592	1 pc	





8104159

Sanibel™ probe tubes for ADI™ ear tips

Fits Corti™, ERO•SCAN™ and the new

OtoRead"	м.		
Part no	Unit		

Box/100 pcs





Probe tip kit for OAE Ø3.3 mm ear tips

Fits all Eclipse™ – use with ADI™ ear tips. Requires no re-calibration.

Part no.	Unit	
8012960	Baa/8 pcs	





Probe tip for Sanibel™ ADI™ ear tips

Fits Titan $^{\text{\tiny{TM}}}$ and easyTymp $^{\text{\tiny{TM}}}$.

Part no.	Unit
8001701	1 pc



Probe Tube Removal Tool

Fits Corti™, ERO•SCAN™ and the new OtoRead™.

Part no.	Unit
8513131	Bag/1 pc





Probe tip kit

Fits AT235™ or newer, AA222™, Sera™,
Tympstar Pro™ Novus™, easyScreen™,
and touchTymp™. Includes probe tip and

Part no.	Unit
8107762	Set/2 pcs





Probe Tube Removal Tool kit

Fits Corti™, ERO•SCAN™ and the new OtoRead™.

Part no.	Unit
8513132	Baa/100 pcs + 1 pc





Sanibel™ probe tip adaptor kit

Converts AT235[™], AT235h[™], AA220[™], AA222[™](blue/new), MI24[™], MI26[™], MI34, MI44[™] for use with ADI[™] ear tips. Includes probe tip, washer, damper. No recalibration required.

Part no.	Unit
8500430	Set/3 pcs





Remote probe tip replacements for ERO•SCAN™

Fits ERO•Scan™.

Part no.	Unit	
8102340	Bag/4 pcs	





Internal probe tip replacements for ERO•SCAN™

Fits ERO•SCAN™.

Part no.	Unit
8102339	Bag/4 pcs



ER10D probe tips for ERO•SCAN™

Pro and OtoRead™

It is recommended to clean the probe tip after every patient and to change the probe tip after 20 tests (40 ears). No calibration necessary.

Part no.	Unit
8012836	Bag/4 pcs





Internal probe tip replacements for ERO•SCAN™ Pro

Fits ERO•Scan™ Pro.

Part no.	Unit
8110364	Bag/4 pcs

Sanibel™ Probe Tube Removal Tool


- A = Probe
- B = Probe Tube Installer
- C = Attachment to probe cable
- D = Tube Remover

To Remove a Probe Tube

- Place the front hole of the Probe
 Tube Tool over the end of the probe,
 as shown. The probe should be
 seated against the face of the tool
 and snap in place.
- 2. Squeeze the tool closed and hold it.
- Twist the tool a couple of times while holding the tool closed and pulling away from the probe.
- 4. The probe tube will pull out from the probe. Discard the probe tube.

To Install a New Probe Tube

- 5. Place a new tube into the hole on the top of the tool.
- 6. Align the probe with the tube.
- Push the probe onto the probe tube until the probe contacts the face of the probe tube tool and a snapping sound is heard.
- 8. Slide the tool off leaving the tube in place in the probe.





3A E-A-RLINK™ foam ear tips

Complies to all insert phones. 13.7mm

Part no.	Description	Unit
8106226	Adult/Medium size	Bag/100 pcs



E-A-RLINK™ infant ear tips

For use with infants or in very small ear canals.

Part no.	Description	Unit
8500300	3.5 mm	Bag/20 pcs
8500310	4.0 mm	Bag/20 pcs



Sanibel™ IP30™ tube connectors For use with IP30™ and E-A-RTone® 3A.

Part no.	Unit
8107391	Baa/10 pcs



3B E-A-RLINK™ foam ear tips

Complies to all insert phones. 9.7mm

Part no.	Description	Unit	
8106227	Pediatric/Small size	Bag/100 p	ocs



Sanibel™ ADI ear tip adaptor

Makes Sanibel™ ADI ear tips compatible with 3A insert phones.

Part no.	Unit	
8500580	Bag/6 pcs	



Tube nipples for E-A-Rtone® 3A

For use with sound tubes 3A.

Part no.	Unit
8500140	Bag/10 pcs



3C E-A-RLINK™ foam ear tips

Complies to all insert phones.

Part no.	Description	Unit
8500100	Jumbo/Large size	Bag/24 pcs
	17.8 mm	



Sound tubes for E-A-Rtone® 3A ABR 50 ohm

Includes adaptor nipples, 1 pc blue and 1 pc red 50 ohm.


Part no.	Unit	
8500720	Set/2 pcs	



Sanibel™ sound tubes and Sanibel™ ADI™ ear tip adaptor kit

Fits E-A-Rtone® 3A inserts and IP30™ insert phones. Includes ADI™ ear tip adaptors, 1 pc blue and 1 pc red 50 ohm.

Part no.	Unit	
8100590	Set/2 pcs	





Sound tubes for E-A-Rtone® 3A 10 ohm

Includes adaptor nipples.

 Part no.
 Unit

 8107396
 Bag/4 pcs




Sanibel™ IP30 and E-A-Rtone® sound tubes and nipples

Fit all IP30 $^{\rm m}$ and E-A-RTone $^{\rm e}$ versions. Includes tube connectors, 1 pc blue and 1 pc red 50 ohm.

 Part no.
 Unit

 8501072
 Set/2 pcs





Sound tubes for E-A-Rtone® 3A 10 ohm

Includes adaptor nipples.

 Part no.
 Unit

 8103605
 Bag/4 pcs





E-A-Rtone® replacement fastener clip

Clip for fastening insert phone on clothes.

 Part no.
 Unit

 8000744
 1 pc



Sanibel™ fastener kit for IP30™

Includes velcro neckstrap and clip.

Part no.	Unit	
8507410	Set	



Sanibel™ ear tip adaptor cleaning tool

Small brush for cleaning the Sanibel[™] probe tip adaptors before disinfection.

 Part no.
 Unit

 8500830
 Bag/10 pcs



Part no.	Description	Unit
8109670	Adult 3A	Kit/100 tests
8109671	Pediatric 3B	Kit/100 tests



Wax Brush 3 in 1

Includes:

- Full brush to remove unwanted wax from hearing aids.
- Magnet at the top of the brush to remove batteries
- Wire look at opposite end

Part no.	Unit
9003569	1 pc



HIF

Temporary Moulds

Used as a temporary/test tip with any Post Auricular Hearing Aid (BTE). Pre-bent Clear #13 tube. Standard 13.7mm diameter.

Part no.	Unit
9002602	Baa/25 pcs







Sanibel™ Infant EarCup™ adaptor set

Connects insert phones with Infant EarCups™. NB: Do not contain filters.

Part no.	Unit
8500450	Set/2 pcs



3M Coban Self-Adherent Wrap

100mm x 2m. Sticks to itself without need for adhesive, pins or clips for fast and easy application.

Part no. Unit
000420 1 Roll



(DPC) DP Probe Couplers

For use with UGD, SND, SGD probes

Part no.	Unit	
DPC-10	Bag/10 pcs	
DPC-100	Bag.100 pcs	





Sanibel™ sound tubes and Sanibel™ Infant EarCup™ adaptor kit

Includes adaptors, 1 blue and 1 red tube.

Part no.	Unit
8100577	Set/2 pcs





Newborn wipes

Fragrance free.

Part no.	Unit
9000287	Canister/ 125 wipes





For use with a range of different providers of newborn hearing screening equipment:

Grason Stadler: Novus. Interacoustics: Titan, Sera. Maico: Easyscreen, MBII Classsic. Madsen/ Otometrics: Accuscreen. Natus: Biologic Abear, Echoscreen.







Use Sanibel™ probe tubes for any equipment

We often hear from the industry that only the manufacturers' specific probe tubes should be used with their specific equipment. However, a recent comparative study by Sanibel Supply® disproves this myth.

In this study, our in-house audiologists and engineers compared the effects of different probe tubes on an Affinity^{2,0} system. They compared AudioScan tubes (red), Otometrics tubes (blue), Sanibel™ tubes (orange) and Sanibel™ tubes (arev).

This involved running each tube via the 'tube calibration' process and then comparing the following REUG curve when the tube had been left in the calibration position.

The investigation showed no significant difference in the REUG curve and all probe tubes achieved the same result when the calibration process had been run correctly.



Silicone tubes for SPL60™

For use with Affinity™ REM probe.

 Part no.
 Unit

 8500700
 Bag/50 pcs



Sanibel™ REM probe tubes, type A

Universal fit, length 76 mm, beveled 45° angle tip and black marker ring.

Part no.	Unit
8104116	Tray/50 pcs



Sanibel™ REM probe tubes, type B

Cost efficient, length 80 mm and marker ring.

Part no.	Unit	<u> </u>
8500360	Bag/40 pcs	_ (



Audioscann probe tubes

Compatible with Audioscan, Axiom and Verifit equipment

Part no.	Unit
8251215	Tray0 pcs







Pre-bent tubing

Part no.	Description	Unit
900020	Standard #13	Bag/10 pcs
900021	Standard #12	Bag/10 pcs
900022	Thick #13	Bag/10 pcs
902085	Extra Thick #13	Bag/10 pcs



Pre-bent tubing

Moisture Free


Part no.	Description	Unit
900037	Standard #13	Bag/10 pcs
900025	Thick #13	Bag/10 pcs
900038	Double Wall #13	Bag/10 pcs



Libby horn tubing

Moisture Free

Part no.	Description	Unit
900032	3mm	Bag/10 pcs
900033	4mm	Baa/10 pcs



Part no.	Description	Unit
900035	Thick #13	Bag/10 pcs
900036	Double Wall #13	Bag/10 pcs



Libby horn tubing

Part no.	Description	Unit
900030	3mm	Bag/10 pcs
900031	4mm	Bag/10 pcs



Stay dry preformed tubing

Part no.	Description	Unit
9002585	Medium	Bag/10 pcs
9002586	Thick	Bag/10 pcs
9002587	Extra Thick	Bag/10 pcs





Brass metal lock

Part no.	Description	Unit
900023	Thick	Bag/10 pcs
900025	13 Thick wall moisture free	Bag/10 pcs







Plastic tube lock with tubing

Plastic lock

Part no.	Description	Unit
900039	13 Thick wall tubing + plastic lock	Bag/10 pcs
900040	Thick Tubing + small plastic lock	Bag/10 pcs
900041	Thick Moisutre Free + plastic lock	Bag/10 pcs
900042	Std Moisture Free + small plastic lock	Bag 10 pcs
900043	Thick-wall + small plastic lock	Bag/10 pcs



Surgical Blades and Handles

Suitable for cutting HA tubing.

Part no.	Description	Unit
8515427	Surgical blade handle No. 3	1 pc
8515428	Surgical blades Noll	100 pcs




Loctite 406 Instant adhesive

An instant adhesive designed especially for the fast bonding of plastics and elastomeric materials where very fast fixturing is required. This product offers exceptionally low viscosity.

Part no.	Unit	
8005067	1 unit	





Tubing Puller

For replacing worn out tubing, Pulls tubes easily through any type of mold or canal tip

Part no.	Unit
9004935	1 pc



Modification Tool Kit

For hard moulds

Part no.	Unit
90005164	1 set



Earlight with Scaled Top

Part no.	Description	
8519357	Luxamed Ear light with	
	acrylic scaled top, black	
8519358	Luxamed Acrylic Top with Scale	



Dremel 200

Comes with 30 different accessories and a 2.4mm Collet.

Part no.	Unit
0006052	1.00+



Modification Tool Kit

For soft moulds

Part no.	Unit
00005165	1 cot



Stethoclip

Hearing Instrument listening stethoscope for reassuring the fitting process

Part no.	Description	Unit
8515600	Stethoclip, black	1 unit
8520346	Ear funnels	1 unit
8520347	ITE Adaptor	l unit



Take care of your hearing aids

Take good care of your hearing aid. Your Hearing instrument is a high tech device needing regular care in order to provide you the best performance every day.



Cleansing tablets

Part no.	Unit
8520148	4 pk



Dehydrating capsules

Part no.	Unit
8520149	20 pk



Drying cup

Part no.	Unit
8520150	1 pc





Cleaning cup

Part no.	Unit
8520151	1 pc



Otofloss

New

Part no.	Unit	
8520152	1 pc	



Perfect Dry LUX box for all types of hearing aids

Part no.	Unit	
8520154	1 pc	



Starter kit

Part no.	Unit		
8520155	kit		



Drying spray

Part no.	Unit	
8520153	1 pc	





Cleaning spray

Part no.	Unit
8520147	1 pc





Protect your ears

Part no.	Description
8520053	Dermoprotective Gel



Sanohra Earplugs Swim

Sanohra Swim earplugs. Enables people to swim without getting water in the ear, reducing ear infections.

Part no.	Description	Unit
8519353	EarPlugs Swim Kids, EN	1 set
8520143	EarPluas Swim Adult, EN	1 set



Cleaning wipes

Part no.	Unit
8520146	30 pk



Protect your ears

Part no.	Description
8520063	Earspray





Protect your ears

Part no.	Description
8520052	Natural lubricant



Sanohra Earplugs Fly

Sanohra Fly ear plugs. Reduces ear pressure experiences when flying.

Part no.	Description	Unit
8519351	EarPlugs Fly Kids, EN	1 set
8519352	EarPlugs Fly Adult, EN	1 set

SANOHRA fly ear plugs

Beyond "noise protection" it can be important to protect the ear drum (tympanic membrane) from changes in external pressure that are too rapid, e.g., during take off and landing phases in airplanes. For airplane passengers experiencing a problem with the pressure equalisation, the use of a special ear protection product means freedom from an unpleasant feeling or piercing pain during flight.

SANOHRA swim ear plugs

In case of injury to the ear drum or insertion of a vent tube, it is necessary to prevent water from reaching the ear drum through the outer auditory canal. Swimming, bathing or showering might cause inflammation of the middle ear. Staying in chlorinated water for an extended period of time can cause the skin in the outer auditory canal to swell and cause inflammation (otitis externa).








Foam Otoblock - Small Pink

Diameter approx, 1.1cm, length, 1.2cm

Part no.	Unit
Otoblkpink	Pkt/100 pcs





Foam Otoblock - Small Grey

Diameter approx, 0.95cm

Part no.	Unit
825403202	Pkt/50 pcs



Cotton Otoblock - Medium White

With white safety thread. Diameter approx, 0.75cm

Part no.	Unit
9140798	Pkt/50 pcs





Foam Otoblock - Medium White

Diameter approx, 1.1cm

Part no.	Unit
Otoblkwhite	Pkt/50 pcs



Foam Otoblock - Standard Yellow

Cone shaped dams with thread. Comes in standard size for CIC Impressions. Diameter approx 1.1 cm, length 0.7cm

Part no.	Unit
Otoblkyellow	Pkt/100 pcs



Cotton Otoblock - Large Red

With red safety thread. Diameter approx, 0.85cm

Part no.	Unit
9140799	Pkt/50 pcs






Foam Otoblock - Large Blue

Diameter approx, 1.25cm

Part no.	Unit
Otoblklarge	Pkt/50 pcs





Cotton Otoblock - Small Blue

With blue safety thread. Diameter approx, 0.65cm

Part no.	Unit	
9140797	Pkt/50 pcs	





Professional Syringes

Part no	Description	Unit	
8900102709	4mm Grey	1 pc	
923000002	3mm Blue	1 pc	



Bite Blocks

Made from styrofoam these diposable bite blocks are used to hold the patient's mouth open in a comfortable position when taking open jaw ear impressions

Part no.	Unit	
8519356	Bag/100 pcs	







Mini Cleaning Brush

B & C components for ear impressions

Part no.	Unit
8251805605	1 nc



Earmold Air Blower

Part no.	Unit	
8525562	1 pc	



OTO-EASE™ Earmold & Hearing **Aid Lubricant**


A grease-less bacteria free liquid lubricant which eases insertion of ear-molds, I.T.E and C.I.C hearing aids

Part no.	Unit
9003259s	1x 15ml bottle



Flat Type Spatula





Bite Blocks

Made from styrofoam these diposable bite blocks are used to hold the patient's mouth open in a comfortable position when taking open jaw ear impressions

Part no.	Unit
8519356	Bag/100 pcs





A-Zoft Ear Impression Material

B & C components for ear impressions

Part no.	Unit
87825015004	2x 800g jars



Dreve 4545 Otoform AK Impression Material Kit

8 x S 50 double cartridges of 2 x 25ml A+B

Part no.	Unit
9004439	2x25ml



A-Zoft Ear Impression Material

B & C components for ear impressions

Part no.	Unit
87825015105	40 foil bags



Otoform A/K 2 part silicone

Proven kneadable impression material, non-sticky and smooth with a slightly pressure-forming consistency. An even and controlled pressure is applied onto the ear tissue. Otoform AK is advantageous for little resilient ear tissues for fitting and sealing of earmolds

Part no.	Unit
9004438	2 pcs



Detax Pink Putty Impression Kit

This impression putty has a soft creamy consistency for easy mixing and removal. The Kit includes two dosing spoons and 500ml of catalyst and 500ml of base material to produce approximately 60 impressions.

Part no.	Unit
A16738	set





Part no.	Description	Unit
8522108	Injection for 48 ml double cartridges	1 pc
8522434	Injection for 50 ml double cartridges	1 pc



Mixing Tips for 50ml double cartridge

Use with 8522433.

Part no.	Description	Unit
8522432	5.4mm	Baa/100pcs



Mixing Pads

Disposable paper pads for mixing impression material. Foam, no slip surface on bottom of pad. Package of two pads with 35 sheets per pad. Dimensions $6" \times 6"$

Part no.	Unit
8522106	100 sheets



Otoform® A flex

Soft, flexible, easy processing. Recommended for CIC impressions, narrow auditory canals or people sensitive to pain. Final hardness approx. 30Shore A.

Note: requires an impression gun and mixing tip.

Part no.	Description	Unit	
8521806	48 ml double cartridges	8 рс	\$160
8522433	50 ml double cartridges	8 pc	\$170



Mixing Tips for 48 ml double cartridge

Use with 8521806.

Part no.	Description	Unit
8522430	5.4mm	Bag/100pcs



Mixing tray

Comes with dosing scale

Part no.	Unit	
8522105	1 pc	







Sanibel™ snap electrodes

Pre-gelled with a very gentle hydrogel, perfect for ABR testing on infants. Size 2.2 x 3 cm. Item no. 8107137 is packed in 20 smaller pouches of 3 electrodes. Intended use areas: Newborn, diagnostic follow-up, ABR, C-Vemp.

Part no.	Unit
8107137	Box/60 pcs (20x3 pcs)







Sanibel™ tab electrodes

Pre-gelled cloth electrode with a very gentle hydrogel, perfect for ABR testing on infants. Size 2.2 x 3 cm. Packed in 20 smaller pouches of 3 electrodes. Intended use areas: Newborn, diagnostic follow-up, ABR, C-Vemp.

Part no.	Unit
8107410	Box/60 pcs (20x3 pcs)



Sanibel™ vinyl tab electrode

For alligator clip cable. Solid gel. Size 2.5 x 2.2 cm. Packaged in convenient re-sealable pouches.

Part no.	Unit
8107426	Bag/100 pcs (10x10 pcs)






E-A-RLINK™ gold foil tip trodes

For ECochG. Use with cable 8011343.

Part no.	Description	Unit
8500370	Adult 13 mm	Bag/20 pcs
8500440	Pediatric 10 mm	Bag/20 pcs





Sanibel Supply* offers various kinds of surface electrodes for ECG as well as EEG and ABR. Surface electrodes vary in size, adhesion, connectivity and frequency (test category).

The frequency plot to the above shows the electrical stimuli sent directly through the system and compared to

Choose the correct electrode for optimal performance

the surface electrode (here ABR). The two lines are similar and demonstrate that the surface electrode can be used for ABR testing.

At Sanibel Supply*, we always test and validate our surface electrodes to ensure optimal performance for your testing!





Vermed Economy line 3-pack solid gel snap Electrode

Solid gel type snap electrode suitable for newborn screening and diagnostic followup. Latex free.

Intended use areas: Newborn, diagnostic follow-up all ages.

Part no.	Unit
8516252	Bag/50 pcs



Disposable electrodes



Vermed Wetgel Snap Electrode

Pre-gelled hydrogel electrode with Ag/ AgCl sensor. Very suitable for newborn screening or diagnostic follow up. Packaged 3x20 pcs. per pouch. Intended use areas: Newborn, Diagnostic follow-up, ABR, C-Vemp.

Part no.	Unit
8516221	Bag/60 pcs (3x20 pcs)







Vermed Tab Electrode

Rectangular tab electrode for newborn screening. 12,7x22,9 mm. Conductive Adhesive Hydrogel. 10x10 pcs. per pouch. Intended use areas: Newborn, O-Vemp.

Part :	no.	Unit
--------	-----	------

8516225 Bag/1000 pcs (100 pcs per bag)



Vermed Hydrogel Snap Electrode

Large soft electrode size making it very suitable for newborn screening and ABR measurements. Cloth material makes it convenient to wear. Size: 22.2. Intended use areas: Newborn, ABR.

Part no.	Unit
8516907	Bag/60 pcs (3x20 pcs)





Vermed Hydrogel Snap Electrode

Smaller than the rectangular version, but otherwise similar. Good soft electrode size making it very suitable for newborn screening and ABR measurements. Cloth material makes it convenient to wear. Intended use areas: Newborn, ABR.

Part no.	Unit
8102022	Baa/60 pes (3v20 pes)







Vermed Tab Electrode

Rectangular tab electrode for newborn screening. 22.9x34 mm. Conductive Adhesive Hydrogel. 10x10 pcs. per pouch Intended use areas: Newborn, O-Vemp.

Part no.	Unit
8516227	Box/100 pcs





Vermed Wet Gel foam

Electrode with an Ag/AgCl sensor.

10x6 pcs per pouch. Diameter 38.1 mm

Part no.	Unit
8108563	Bag/60 pcs







Vermed Prewired neonatal ecq

electrode for use in pediatric and neonatal settings. 3 electrodes per pouch. square shape 22.2x22.2 mm

Part no.	Unit
8107425	Bag/300 pcs (3x100 set)







Vermed Economy line Snap Wet **Gel Electrode**

Non-woven cloth electrode with conductive wet gel. Passed ISO 10993 and is Latex free.

Intended use areas: Newborn, Diagnostic follow-up all ages.

Part no.	Unit
8516251	Box/50 pcs (20x50 pcs)







Vermed Allround ECG electrode

Provides very good impedance and it well suited for both newborns and adults. Can be used both as snap and tab electrode with alligator cables. Pouches come in 3 pcs, value packs available in 600 pcs. per order. Please contact for order.

Part no.	Unit
8516906	Bag/5 pcs





Disposable electrodes



Ambu® 0415M tab electrode

Medium size tab electrode with solid gel for 12-lead ECG. Presented conveniently with 10 pieces to a card. Fish tail tab design reduces electrode lift-offs. Suitable for both audults and children.

Intended use areas: ABR Diagnostics of all kinds with Tab electrodes.

Part no.	Unit
8516138	Bag/100 pcs (10x10 pcs)





Ambu° Ground electrode 71410-M/1

Part no.	Unit
8500680	1 pc



Ambu[®] Neuroline 720 - with 100 cm cable and k-connector

Ambu Neuroline 720 electrodes are pre-gelled and self adhesive. They are ideal for use in sleep registration where the strong, yet skin-friendly, adhesion ensures that the electrode stays in place for the entire examination, even if the patient sweats. Intended use areas: ABR Diagnostics of all kinds.

Part no.	Unit	
8503244	Bag/10 pcs	



Ambu° 2300 Blue Sensor Tab electrode

The medium size Ambu BlueSensor 2300 electrode features a skin-friendly solid gel that provides optimal signal quality recordings during 12-lead diagnostic resting ECG application. Very suitable for ABR diagnostics of all kinds. Intended use areas: ABR Diagnostics of all kinds

Part no.	Unit
8516139	Box/1000 pcs (10x100 pcs)



Ambu® Neuroline 720

Ambu Neuroline 720 electrodes are pregelled and self adhesive. They are ideal for use in sleep registration where the strong, yet skin-friendly, adhesion ensures that the electrode stays in place for the entire examination, even if the patient sweats. Intended use areas: ABR Diagnostics of all kinds.

Part no.	Unit
8500060	Baa/25 ncs



Ambu[®] Blue Sensor™ N-10-A

Foam ECG electrode with 10 cm wire and 4 mm connector, for use with universal electrode cable. Wetgel.

Part no.	Unit
8500670	Bag/25 pcs



Ambu[®] Blue Sensor™ N-00-S

Foam ECG electrode for use with button cables. Wetgel.

Part no.	Unit	
8013029	Bag/25 pcs	



Ambu° Blue Sensor™ P-00-S

High quality Ag/AgCl sensor – ensures excellent trace quality.

Part no.	Unit
8504783	Bag/50 pcs



Ambu[®] Neurology Needle Electrode

This needle electrode features exceedingly low penetration-resistance levels. The needles have been coated with smooth coating which reduces friction when inserting and working in the target area.

Part no.	Unit	
8519037	Box/40 pcs	


Disposable electrodes



Ambu® Blue Sensor™ M-00-S

High conductive wet gel with occlusive backing material making it liquid resistant

Part no.	Unit	
9000801	Baa/50 pcs	



Ambu[®] Blue Sensor™ N-50-K

High quality Ag/AgCl sensor – ensures excellent trace quality.

Part no.	Unit
8521842	Bag/25 pcs





Ambu® Blue Sensor™ NF

Specifically developed for babies and smaller children.

Part no.	Unit	
9018110	Baa/12 pcs	





Ambu White Sesnsor 40554

Specifically designed for babies and newborns.

Part no.	Unit
8520431	Bag/3 pcs



Surface electrodes tailored to your need

At Diatec we believe that it is important to choose the right surface electrode for each specific recording and patient type. It is important to use electrodes with similar chloride salt content when testing.

Wetgel or hydrogel electrodes?

A wetgel electrode consists of water, thickener, ionic salts, surfactants and bactericides or fungicides. Wetgel migrates easily and rapidly onto the skin surface, reducing impedance levels, which in turn will generally result in lower noise in the recording.

A hydrogel surface electrode is less adhesive when compared to wetgel, making it particulary suitable for newborns.

Chloride salt concentration

The higher the chloride salt content, the more conductive the electrode. Higher salt content (10% and above), pregelled surface electrodes are useful for making fast, high quality measurements of biopotentials. High conductivity electrodes generally have reduced artifact, due to the lower impedance generated between electrode and skin surface. As the chloride salt content drops, there is a corresponding reduction in electrode conductivity.

If the chloride content is below 10%, the electrode can then be employed for longer recording times, with a reduced risk of skin irritation. As a rule of thumb, wetgel has a chloride content of 10% and hydrogel 4%.

The impedance of the electrode/skin junction is highly dependent on the electrolyte type and the chloride salt concentration. For example, a hydrogel electrode with 4% chloride concentration will have approximately 10x higher impedance than a wet liquid gel electrode with 10% chloride concentration, after first application to the skin.

Adhesion

An advantage of a conductive adhesive is that the active surface can be as large as the electrode itself because no additional adhesive strip is needed. The Sanibel disposable electrodes come with a range of adhesive qualities. All electrodes are designed to adhere well to skin surfaces, but the adhesion characteristics can be grouped according to the application:

- Strong adhesive
- Moderate adhesive

Strong adhesion is not recommended for newborns.

Skin preparation

For optimum electrode to skin conductivity, the skin should be lightly brushed with a gentle abrasive gel/wipe.

An alcohol wipe is not recommended to improve conductivity as this will only dry out the skin surface.





Sanibel™ ECochG starter kit

For Eclipse $^{\rm T\!M}$ and Bio-logic $^{\rm S\!M}$ AEP system. Includes cable, gel and 2 pcs TM electrodes.

Part no.	Unit	
8013209	Set/4 pcs	



Sanibel™ ECochG cable

For Eclipse™ and Bio-logic® AEP system. 1.5 mm/.060″ female connector. For use with Sanibel™ TM electrodes.

Part no.	Unit	
8500650	1 pc	

After use ...

After use, cables should be wiped with a soft cloth that has been dipped in a mild detergent solution or alcohol.

Contact between the cleaning solution and the connectors should be avoided.



Disc electrodes

For use with Ten20™ electrode paste. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable.

Part no.	Unit
8011498	Set/3 pcs RWB
8011345	Set/5 pcs RWBYK



Gold-plated cup electrodes

Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). Ø 10 mm gold-plated. 50 cm cable.

Part no.	Unit
8107418	Set/3 pcs RWB
8107145	Set/5 pcs RWBYK



ABR

Alligator clip cable

For use with tab electrodes. 1.5mm/.060" female connector. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K).

Part no.	Unit
8011412	Set/3 pcs RWB
8500900	Set/3 pcs WYK
8011352	Set/5 pcs RWBYK



For use with TTE25 gold foil tip trodes. 1,5 mm with touch proof connectors and plated clamp clips. 50 cm long.

Part no.	Unit
8011343	Set/2 pcs



Universal electrode cable

Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable.

Non-tangling.

Part no.	Unit	
8011353	Set/3 pcs RWB	
8011344	Set/5 pcs RWBYK	



For use with snap electrodes. Available in Red (R), White (W), Blue (B), Yellow (Y),

Black (K). 50 cm cable. Non-tangling.

Part no.	Unit
8011354	Set/3 pcs RWB
8011346	Set/5 pcs RWBYK



Alligator clip leads

For use with tab style disposable electrodes

Part no.	Description	Unit
301633	R/WT/BL, 18"	Set/3 pcs
9417000	R/G/BK, 1.25m	Set/3 pcs
301635	R/WT/BL/G/BK, 18"	Set/5 pcs



Pinch style cable

For use with snap electrodes. 1.5mm/ .060" female connector. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K).

Part no.	Unit
8500390	Set/3 pcs RWB
8519048	Set/3 pcs WYK
8500410	Set/5 pcs RWBYK



ABR

Sanibel™ EP electrode jumper cable

Allows two of the preamplifier inputs to be connected, e.g. for linking mastoid electrodes for cortical testing.

Part no.	Unit	
8107872	og [



ABR

Pinch clip leads

Lead wires for use with snap disposable electrodes

Part no.	Description	Unit
301663	R/WT/BL, 18"	Set/3 pcs
301665	R/WT/BL/G/BK/, 18"	Set/5 pcs
301666	R/WT/BL/G/BK/Y, 18"	Set/6 pcs





NuPrep™ skin prep gel, 114 g

Efficient, yet gentle skin abrasive. Chloride and acetone free.

Part no.	Unit	
8500320	Box/3 tubes	



Lectron® II electrode gel, 250 ml

Highly conductive, chloride-free electrode gel. Recommended for use with Sanibel™ TM electrodes for ECochG.

Part no.	Unit	
8102626	Tube/250 ml	



Sonavelle® electrode gel, 250 ml

Formaldehyde-free and glycerol-free.

Part no.	Unit
8500380	Bottle/250 ml



Spectra 360[®] electrode gel, 250 g

Salt- & chloride-free conductive electrode gel suitable for long-term applications.

Part no.	Unit
8500740	Tube/250 a



Everi™ abrasive paste, 160 g

Reduces the skin impedance during EEG, EP, ECG and Holter Stress test examination.

Part no.	Unit
8012952	Bottle/160 g



Ten20™ conductive EEG paste, 114 g/228 g

Adhesive gel for disc electrodes. Available in pack of 3x4 oz/114 g tubes, 3x8 oz/228 g jar.

Part no.	Unit
8012924	Box/3 tubes
8500330	Box/3 jars



SignaGel®, 250 g

Hypoallergenic, bacteriostatic conductive saline electrode gel.

Part no.	Unit	
8102628	Tube/250 g	



Niltac™ sting-free adhesive remover

The perfect solution to removing tape and bandages without removing hair or damaging sensitive skin.

Part no.	Unit	
91706525	150ml	



3M™ Red Dot™ trace prep, 5 m

One-step skin preparation. 18 mm x 5 m / 0.70" x 196".

Part no.	Unit
8500030	Dispenser/5 m



3M™ Transpore™ Surgical tape, 9.14 m

Transparent, latex-free, water-resistant and porous for breathability. 25 mm x 9.14 m / 1" x 10 LY.

 Part no.
 Unit

 8500750
 Box/12 rolls



Gauze swabs

Non-sterile, 7.5 x 7.5 cm and 8 layered.

 Part no.
 Unit

 8500040
 Pack/100 pcs.



Small Alcohol Pads

70% isopropyl alcohol. Breimer swabs 30 x 110mm

Part no.	Unit	
9581200	Box/200 pcs	



3M™ Micropore™ Surgical tape, 9.14 m

Latex-free and hypoallergenic. Holds well on damp skin. 25 mm x 9.14 m / 1" x 10 LY.

Part no.	Unit	
8500760	Box/12 rolls	



Q-tips

Cotton tipped applicators for patient preparation. Non-sterile. Each pack contains 100 pcs.

Part no.	Unit	
8500810	Box/10 packs	





Uni-Solve

Adhesive remover wipes

Part no.	Unit
90402300	Box/50 wipes





3" Disposable ear cushion cover

Replaces disinfactants. Headphone Covers for classic audiometric headphones.

Part no.	Unit
8107419	Box/500 pcs









Probe cleaning tool

For Interacoustics Titan™/AT235™/AA222™/Sera™/MT10™/Amplivox Otowaves/MAICO Mi XX + touchTymp $^{\text{TM}}$ and GSI TympStar $^{\text{TM}}$ Pro.

Part no.	Unit	
8011351	1 pc	



Disinfectant detergent, 1000 ml

Alcohol-free, ready to use solution for a quick and extensive disinfection of non-invasive medical devices.

Part no.	Unit
8501706	Bottle, 1000 ml



4.5" Disposable ear cushion cover

Replaces disinfactants. Headphone Covers for noise exclosure audiometric headphones. For HDA200/DD450/DD65v2/Peltor headsets.

Part no.	Unit	
8513404	Baa/100 pcs	



Keyboard protective foil

Place a protective foil over your keyboard to minimize contamination with grime, bacteria and germs, or liquid.

Part no.	Unit	
8500460	1 pc	





Bridge & Implant floss

Used to clean probe tips.

Part no. Unit

Box/30 pcs 8003566



3-in-1 floss

Used to clean probe tips.

Part no. Unit

Box/100 pcs 8003565



ABR IMP OAE
Spiraflex Cleaners

Used to clean probe tips.

Part no. Unit

8902124100 Bag/30 pcs



Large alcohol swabs

70% isopropyl alcohol, 2% CHG. PDI swabs. 150 x 162 mm

Part no. Unit

900331 Box/200 pcs





PDI Alcohol wipes canister

70% isopropyl alcohol

Part no.	Unit
900325	200 sheets





Spirigel sanitising hand gel, pump

Kills germs without water, air dries in seconds, added moisturiser.

Part no.	Unit	
9045940	500ml	



PDI Detergent wipes, bucket

Wipe size: 300 x 245mm

Part no.	Unit	
9000175	Bucket/225 pcs	



PDI Detergent wipes, refill

Solo lid bag. Wipe size: 300 x 245mm.

Part no.	Unit
8000210	Bag/225 pcs



Infection **control**

A clean clinic is also part of your professional image

Infection control is so much more than cleaning probe tips, protecting your personal computer's keyboard, and maintaining high standards within hygiene using e.g. disinfectant detergent for a quick and extensive disinfection of non-invasive medical devices.

When it comes to infection control, preventing germs, dirt and dust from spreading via other sources, is also of critical importance – both to the safety of your patients and yourself! Gloves, shoe covers, facemasks, and aprons are just some of the "accessories", which can help you succeed in this.

Maintain high, professional standards and make your clinic an ever-safe environment for your patients, making the wearing of gloves, facemasks and other body protection wear a natural part of the "patient interaction" process.

Personal Protective Equipment





Bouffant Caps

Perform high Hygiene practices in front of your patients or clients.

Part no.	Unit
9011312	Bag/250 pcs (21 Inches)
9250112	Bag/250 pcs (24 Inches)



Full Face Shield

Anti-fog, Anti-glare, Anti-static

7 11 10 10 9, 7	and grand, rand
Part no.	Unit
9001063	1 pc





Clear safety glasses

Part no.	Unit	
9000564	1 pc	



Disposable facemasks with loops

Perform high Hygiene practices in front of your patients or clients.

Part no.	Unit
991003	Pack/50 pcs



Disposable apron

Keep a professional appearance also after servicing a hearing instrument.

Part no.	Unit
8519348	Box/100 pcs





Promed soft examination gloves

Latex and powder-free and provide excellent protection. Our gloves are available in many sizes and will accommodate most hands

Part no.	Description	Unit
9173070	Small	Box/200 pcs
9173071	Medium	Box/200 pcs
9173072	Large	Box/200 pcs
9173073	Extra Large	Box/200 pcs



Disposable Shoe Covers

Perform high Hygiene practices in front of your patients or clients.

Part no.	Unit
8519349	Box/100 pcs (Regular)





Disposable VNG goggle foam pads

For Interacoustics $^{\! *}$ VNG mask.

 Part no.
 Unit


 8011527
 Box/24 pcs





Disposable specula for Difra Airstar

Part no.	Description	Unit
8011457	2.5 mm spec. (green)	bag/100pcs
8011456	4 mm spec (blue)	bag/100pcs




BAL

Water irrigator tips

For use with $\mathsf{Interacoustics}^*$ Aqua $\mathsf{Stim}^{\mathsf{TM}}.$

Part no.	Unit
8013018	Bag/50 pcs







Disposable ear specula

For AirFx™.







Accessories for Otoscopy

Stylish wall dispenser for ear funnels. Can contain 6 different sizes.

Part no.	Description
8515678	Specula wall dispenser



Video Otoscopes specula

Fits the new Interacoustics* Viot and OriVision otoscope.

Part no.	Description	Unit
8515083	KaWe 2.5 mm speculum	Bag/100 pcs
8511728	KaWe 4 mm speculum	Bag/100 pcs





Disposable speculas for video otoscope

Fits the AURICAL $\!^{\rm IM}$ OTOcam 300, Interacoustics $\!^*$ VIOT v1 and XION $\!^{\rm IM}$ video otoscopes.

Part no.	Unit
8500470	Bag/50 pcs





BAL

Disposable ear specula

For Welch Allyn.

Part no.	Description	Unit
8107750	2.5mm	Bag/250 pcs
8108252	4.0mm	Bag/250 pcs



Disposable ear specula

For Heine Otoscope. Suitable for adult and children.

Part no.	Description	Unit
8450102907	2.5 mm	Bag/50 pcs
8450102806	4.0mm	Bag/50 pcs



Bulb for Heine Mini Otoscope

2.5v Halogen replacement bulb 110

Part no.	Unit	
8900204705	1 pc	

Heine Mini Otoscope



Heine Mini Otoscope

Heine Mini 3000 set. Includes 5 speculas.

Part no. Unit 8900400105 Box Set



Welch Allyn Diagnostic Pocket LED Otoscope

Part no.	Description	Uni	
9000317	Vanilla	1 pc	
9000318	Onyx	1 pc	
9000319	Plum	1 pc	
9000320	Blueberry	1 pc	

Features

- Light-weight Welch Allyn otoscope for comfort and control during examinations
- Portable, to help you meet the challenges you face when treating patients
- SureColour™ LED lamps come standard
- Fibre-optics provide cool light with no reflections or obstructions



Luxamed Otoscope Auris LED 2.5V

Otoscope from German Luxamed. LED light to ensure longitivity and great light source, dimable. Fits Heine/KaWe compatible speculums.

Part no.	Unit
8519433	1 nc

Luxamed Otoscope Auris LED 3.7V USB/PC Charged

Otoscope from German Luxamed. LED light to ensure longitivity and great light source. USB charging. Fits Heine/KaWe compatible speculums.

Part no.	Unit	
8519434	1 pc	
	NEW	5 pcs



Cartoon faces that are placed on the otoscope. Make ear examinations less scary for kids.

Part no.	Unit
8011501	Baa/5 pcs (assorted)









Round nose scissors

Part no.	Unit	
9100234	1 pc	



Tube expander

Part no.	Unit	
930421	1 pc	



Cleaning tool

Part no.	Unit
8902126000	Bag/20 pcs



Nelson tool

Part no.	Unit		
00043331	Sot		



Mini pliers

Part no.	Unit	
9100235	1 pc	



Mini wire cutters

Part no.	Unit	
9100236	1 pc	



Tweezers

Part no.	Unit
9100237	1pc



Tweezer Set

Part no.	Unit		
00005300	Cot		

Paper



Thermal paper for ERO•SCAN™/ OtoRead™, 25 m

Width: 57 mm

Inner diameter: 12 mm Outer diameter: 45 mm

Part no. Unit 8500200 1 roll/25 m





Thermal paper for Sanibel™ MPT-II printer, 20 m

For the Sanibel™ MPT-II printer which is suplied with MAICO® ERO•SCAN™, Grason-Stadler® Corti™ and Interacoustics® OtoRead™.

Width: 57 mm

Inner diameter: 9 mm Outer diameter: 40 mm

Unit Part no.

1 roll/20 m 8029305







AP1300 thermal paper for Titan™, 10 m

Width: 57 mm

Inner diameter: 12 mm Outer diameter: 35 mm

Unit Part no.

1 roll/10 m 8500230









TPR10™ thermal paper for MTP10™,

Inner diameter: 12 mm Outer diameter: 45 mm

Part no. Unit

8500220 1 roll/25 m





Sanibel™ MPT-III thermal paper, 20 m

For Sanibel™ MPT-III printer which is supplied with Interacoustics® AT235™.

Width: 79 mm

Inner diameter: 12 mm Outer diameter: 45 mm

Part no. Unit 8500610 1 roll/20 m





TPR26™ thermal paper for Interacoustics®, 25 m

Fits most popular audiometers and tympanometers.

Width: 112 mm


Inner diameter: 12 mm Outer diameter: 45 mm

Part no. Unit 8500210 1 roll/25 m





Paper



Adhesive thermal paper LPQ80, 195 labels

For GSI TympStar™ Pro. Width: 82 mm Inner diameter: 25 mm Outer diameter: 88 mm

Part no. Unit

8509894 1 roll/195 labels





For OtoRead™ PRT MPT-II printer, 10 m

Part no.	Unit	
8013009	1 roll/10 m	IMP



GSI Thermal paper

For GSI39

Width: 111 x 51 x 12 mm

Part no.	Unit	
8012939	1 roll	



Thermal paper LPQ80, 80 m

For GSI TympStar™ Pro. Width: 80 mm Inner diameter: 13 mm Outer diameter: 80 mm

Part no. Unit 8509893 1 roll/80 m



Adhesive Label Paper roll -HM-E200

Part no.	Unit	
8519836	1 roll/40 labels	 IN



GSI Thermal paper with label

For GSI37

Width: 58 mm

Part no.	Unit	
8012950	1 roll	



Thermal paper for GSI70

For GSI 70.

Width: 109.5 mm

Part no.	Unit
8500240	1 roll



Thermal paper for GSI

117 Thick, Adhesive, 9m long. For GSI Tympstar

Part no.	Unit	
8012937	1 roll	




Audioscan printer paper

For VFI, Axiom and RM500SL Width: 80 mm

Part no.	Unit	
8521224	1 roll	



Wax removal accessories



Safe Ear Curettes allows safe and gently removal of different types of cerumen. Single use ear currettes. Ensure your patients' safety by using these affordable ear currettes. Available as variety pack or single size items. Safe Ear Curettes allow healthcare provides to safely and gently remove cerumen with styles to match different types of wax and personal preference.

Part no.	Description	Unit
8515585	Green MicroLoop®	50 pcs/box
8515586	Orange ControLoop®	50 pcs/box
8515587	Red AngleLoop*	50 pcs/box
8519359	Purple VersaLoop*	50 pcs/box
8519360	Blue InfantScoop®	50 pcs/box
8519372	Yellow CeraSpoon*	50 pcs/box
8519361	White FlexLoop®	50 pcs/box
8519362	Variety Kit, (10 each style	e & 15 white)



Lighted Foreceps


Bionix Lighted Forceps combine illumination and magnification to help healthcare providers take less time to see and safely remove ear canal and nasal passage obstructions.

Part no.	Description
8519392	10 foreceps, 1 light and 1 magnifying lens



Lighted Ear Curette Variety Pack

Part no.	Description
8519363	with 1 Light source and magnifying lens includes all tip
	styles except CeraPik™



Bionix Lighted Ear Curettes

Box of 50 pcs. Use the Bionix lighted ear currettes for safe removal of earwax. The Bionix currettes comes in a variety of different shapes for you to option for your optimal choice in a situation. Bionix Lighted Ear Curettes combine illumination and magnification to help healthcare providers take less time to see and safely remove cerumen from the ear canal.

All versions including 1 Light Source and 1 magnifying lens.

Part no.	Description
8519364	Lighted AngleLoop*, 4mm diameter
8519365	Lighted MicroLoop®, 3mm diameter
8519366	Lighted WaveCurette™, 2.5mm wide tip
8519367	Lighted FlexLoop®, 4mm wide tip
8519368	Lighted CeraSpoon®, 4mm wide tip
8519369	Lighted InfantScoop*, 2mm wide tip
8519373	Lighted VersaLoop®, 3mm wide tip
8519370	Lighted CeraPik™, 2.5mm concave tip. Great for hard, impacted cerumen and foreign bodies.
	8519364 8519365 8519366 8519367 8519368 8519369 8519373



Lighted articulating curettes

The lighted articulated curette combines illumination, magnification and a 90° articulating tip.

Part no.	Unit
9008111	Box/50 pcs

Wax removal accessories



Lighted suction handles

The lighted suction combines illumination and magnification with the efficiency of a suction device. Comes with 25 handles, 25 rubber ends, 1 light source, 1 lens.

Part no.	Unit
9002625	Box/25 pcs



Hal Hen LED Earlight

Part no.	Unit		
9125718	1 set		



EarWay Pro - Wax Removal Helix


Very easy and safe way to remove ear wax. The helix grabs around the ear wax lumb and pulls it out in a comfortable way.

Part no.	Description	Unit
8518025	Pro M, 5.4 mm, Purple,	25 units
8518024	Pro S. 4.2 mm, Blue,	25 units



Earlight with Scaled Top

Part no.	Description
8519357	Luxamed Ear light black, acrylic scaled top
8519358	Luxamed Acrylic Top with Scale







Earlight for wax removal

Part no.	Description	Unit
900106	Set	1 set
900107	STD Straight Tip	Bag/4 pcs
900108	Pros Tip	Bag/4 pcs
900109	Spoon Tip	Bag/4 pcs
900110	Stainless Steel Loop	Bag/2 pcs



Otoclear Irrigation Tips

Bionix OtoClear Ear Irrigation helps provide safe and comfortable ear irrigation without the backsplash, at half the cost of a self-constructed device.

Part no.	Unit	
8519375	Box 40 pcs	



Clean ears spray, 30ml

Dissolves earwax and irrigates ear canal. Simple easy-to-use spray which is suitable for all ages.

Part no.	Unit	
9001258	1 pc	



Otoclear Spray wash Kit

OtoClear Ear Irrigation system directs water to the ear canal wall gently with three divergent streams, instead of directly at the Tympanic Membrane

Part no.	Description
8519371	box 20 irrigator tips, 1 470ml spray wash bottle, 1x ear basin



Otoclear® ear basin

Part no.	Unit
8519390	Box/3 pcs





New

Minston Headlight and Loupe

Minstons Headlight and Loupe combine to give an ultra lightweight and comfortable diagnostic experience. The adjustable magnification of 2.5 - 3.5 x makes it suitable for a wide range of applications including ENT, general examination and minor procedures

Part no. Description

9006020 Minston Headlight

Features

- d Ergonomic Design
- d Variable magnification 2.5-3.5x
- d 10,000-40,000 Lux
- d Patented Heat Dispensing Technology
- d L Shaped Lamphead design
- d Ultra Light & Free
- d Even-light
- d Brightness control
- d Dual Battery
- d Easily adjustable





Micro Suction Products

Otopront VAC18 is a transportable suction unit, specific designed for ear wax removal or other temporary suction purposes within hospitals and clinics.

Part no.	Description
8515194	Otopront VAC 18 suction unit

Accessories for Otopront Suction units

To optimize the use of the Otopront VAC18. Please use the appliances below per demand.

Part no.	Description	Unit
8517966	Suction container system	1 pc
(to be ad	Ided the VAC 18)	
8515372	Refill bags for disposables suction system	1 pc
8520345	Silicon suction hose medical grade 50m	1 pc
8515195	Air filter VAC 18	1 pc
8515197	Internal filter VAC 18 / VAC 28	1 pc
8515199	Lid for suction container VAC 18 / VAC 28	1 pc
8522023	Turmakin Ear Speculum 5.5mm	1 pc
8522024	Turmakin Ear Speculum 7.0mm	1 pc
8522025	Suction cannula 2mm x 80 mm	1 pc















Tongye MC600A Suction Unit 240V

The MC600A Suction Unit is a high quality AC powered suction unit featuring step less suction adjustment, overflow protection and anti-bacterial filter. This unit is suitable for many applications such as aural micro suction, resuscitation and home care.

Part no.	Description
900600	MC600A Suction Unit



Tongye Disposable Suction Unit Liners

To suit MC600 units. Provides a disposable and cleaning free solution for the collection of many different matter and is situated inside the hard reusable cannister included with the Tongye MC600 series suction units.

Part no.	Description
90006001	Disposable suction unit liners



Tongye Overflow Stopper

Overflow stopper for Tongye Suction Units to fit all MC600 Tongye suction units

Part no.	Description
90006003	Overflow stopper



Suction Tube 3M

A 3 meter silicone tube to suit the MC600 series suction units.

Part no.	Description
90006002	5mm x 3m



Tongye Reusable Suction Canniste

Provides collection of many different matter and is suitable for use with the MC600 series suction units

Part no.	Description
90006004	Suction cannister





Tongye Hose Barb

Part no. Description
90006005 Hose Barb for MC600



Suction Controller Sterile

A small piece that goes in line between a suction tube and silicone tube that allows contol of suction pressure at the fingertip

Part no.	Unit
90006006	Baa/10 pcs



Tongye Filter

Part no.	Description
90006007	Filter for MC600



Ear Speculum

Speculums have an Oval shaped tip to match the shape of the ear canal, with a thin 0.5mm wall thickness. Light grey in colour which matches the inside stainless steel speculum and is the best balance between light reflection and contrast.


Part no.	Description	Unit
910004	Extra Small, diameter 3 x 3.5mm	Bag/100 pcs
910003	Small, diameter 4 x 4.5mm	Bag/100 pcs
910002	Medium, diameter 5 x 5.5mm	Bag/100 pcs
910001	Large, diameter 6 x 6.5mm	Bag/100 pcs
910000	Extra Large, diameter 7.5 x 8mm	Bag/100 pcs



Ear Speculum Gruber

Speculums have an Oval shaped tip to match the shape of the ear canal. Chrome plated

Part no.	Diameter - Height	Unit
910005	25mm - 36mm	Bag/4 pcs
	26mm - 39mm	0. 1
	25mm - 39mm	
	26mm - 38mm	



Disposable Suction Tube

A single use disposable tube for the suction of fluid and other media from the ear

Part no.	Description	Unit
90006008	14 Gauge Bent, non-sterile Outside diameter 2.1mm	Bag/100 pcs
90006009	16 Gauge Bent, non-sterile Outside diameter 1.6mm	Bag/100 pcs
90006010	18 Gauge Bent, non-sterile Outside diameter 1.2mm	Bag/100 pcs
9006011	20 Gauge Bent, non-sterile Outside diameter 0.9mm	Bag/100 pcs
9006012	22 Gauge Bent, non-sterile Outside diameter 0.7mm	Bag/100 pcs
9006013	24 Gauge Bent, non-sterile Outside diameter 0.5mm	Bag/100 pcs

Parts for audiometry





HDA200 replacement cable, symmetric

For HDA200 and DD65 headphones. 2 m cable and 6.3 mm/ $30^{\rm o}$

jaak plag.		
Part no.	Unit	
8004172	1 pc	



B71 Bone conductor replacement cable

2 m cable and 6.3 mm/ 30° jack plug.

Part no.	Unit	
8508473	1 pc	



Clothing clip kit

Part no.	Unit
3103992	1 pc



TDH39™/DD45™ replacement cable, symmetric

For TDH39™ and DD45™. 30° jack plug.

Part no.	Description	Unit	
8507203	2 m	1 pc	
8507204	5 m	1 pc	



TDH39[™]/DD45[™] contra headset replacement cable

For TDH39™ and DD45™. Single mono. 30° jack plug.

Part no.	Unit	
8010826	1 pc	







IP30 Contra Insert earphone cable

Cable length: 2 m. Plug type: 30°, 1x6.3 mm gold-plated mono jack.

Part no.	Unit	
8107358	1 pc	







RadioEar Foldable audiometric headband

Black, foldable.

Part no.	Unit	
8010882	1 set	AUD



RadioEar RE-7 Headband

Black.

Part no.	Unit	
8517379	1 set	AUD



Peltor® headset/earmuff replacement hygiene kit

For DD450, HDA200 and DD65.

Part no.	Unit		
8500270	1 set	AUD	



RadioEar DD45 Audiometric headset

Headband: HB7. Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

o. Unit	
54 1 set	





RadioEar DD45 Adult audiometric headset

Headband: HB3045 .Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

Part no.	Unit		
8105268	1 set	AUD	







Part no. Unit 8510117 1 set



RadioEar DD45 transducer

Audiometric earphones engineered to match the performance and function of the most popular earphones available for audiological testing. 10 ohm receiver. Easy installation. Reliable. RoHS compatible.

Part no.	Unit	
8507268	l set	AUD



DD 51 cushion

For DD45™ and TDH39™ headsets.

Part no.	Unit
8010894	1 pc



cushion red & blue, model 51

Medical grade

For TDH39, 49, 50 & DD45 earphones.

101 101100	, 40, 00 Q L	D40 carpilone
Part no.	Unit	
8107611	1 Set	







Foam pad for contra headset

Compatible with e.g. TDH39c $^{\text{\tiny{TM}}}$ and DD45c $^{\text{\tiny{TM}}}.$

Part no.	Unit	
8507292	Bag/10 pcs	AUD



RadioEar replacement pads (3288-1)

For RadioEar Spring Steel headband.

Part no.	Unit			
8504777	Bag/10 pcs	AUD	HIF	ABR



Ear cushion HZP 09 pair

Ear cushion foam (1 pair) suitable for: PC 20, PC 25, PC 30, PC 31, PC 35, PC 36, PC 130, PC 131, PC 135, PC 136, PC 140, PC 141, PC 145, PC 146.

Part no.	Unit	
8104416	1 pair	AUD



Headband cushion

For RadioEar HBA foldable audiometric headband.

Part no.	Unit	
8010886	1 pc	AUD



RadioEar B71W Bone transducer headset

Ohm: 10 ohm. Cable length: 79"/200 cm molded cord. Plug type: 30°, lx6.3 mm goldplated mono jack. Headband: P3333

Part no.	Unit		
8510370	1 pc	AUD	



RadioEar Spring Steel headband (P3333)

Part	no.	Unit			
8507	'213	1 pc	AUD	HIF	ABR



Headband cushion

For Sennheiser HDA200, RadioEar DD450, DD65v2 and Peltor audiometric headbands.

Part no.	Unit	
8000770	1 pc	AUD



RadioEar AMBAND bone conductor headband (P1188)

For adult, children and infants.

Part no.	Unit			
3507214	Bag/10 pcs	AUD	HIF	ABF



Replacement cable for bone conductor B71/B81 (with safety fixation point)





OAE Probe body nozzle

Part no.	Unit
580-PRBNZL	1 pc



Ultra-Vac Filter (old)

Part no.	Unit	
8111312	1 pc	



Ultra-Vac 4 piece tool kit

Part no.	Unit	
8111307	set	



Ultra-Vac Needle pack

Part no.	Unit
8111306	Bag/6 pcs



Ultra-Vac Filter (new)

Part no.	Unit	
8506932	1 pc	





Ear coupler sealing wax

Use to create a seal between coupler and hearing instrument.

Part no.	Unit	
8500080	1 pc	





Battery tester

Part no.	Description	Unit
8515452	Battery Tester	1 pc
8520138	Battery "carrier"	1 pc



Lithium battery 3.7V

For Interacoustics® Callisto™ and Titan™.

Part no.	Unit		
8104365	1 pc		







Welch Allyn 3.5V Lamp

For Welch Allyn Macro-View otoscope.

Part no.	Unit
940065000	1 pc



Battery Adaptor


For Interacoustics Affinity, $^{\text{TM}}$ Callisto $^{\text{TMM}}$, Equino x^{TM} , MedRx Avant Hit+ and Siemens Unity 2.

Part no.	Description	Unit
8011198	BAA675	1 pc
8011199	BAA13	1 pc
8011200	BAA312	1 pc
8011201	BAA10	1 pc
8011202	BAA5	1 pc



To order, please send us an e-mail info@diatec-diagnostics.co.nz or contact us on 09 279 8833

diatec



Diatec New Zealand• info@diatec-diagnostics.co.nz Diatec is a company within the William Demant Group

Disclaimer

This catalogue has been published for worldwide distribution. Some parts described may NOT BE AVAILABLE in your area. Your local distributor will be pleased to give you details.

For advice on compatibility and instructions of use, please consult your operation manual or local distributor. Sanibel Supply A/S assumes no legal liability or responsibility for the accuracy, typing errors or mistakes, correctness or completeness of any information in this flyer. The information in this flyer was correct to the best of our knowledge at the time of printing. Sanibel Supply A/S reserves the right to make changes at any time in specification of products described in this flyer without notice or obligation.

MADSEN Accusoreen' and MADSEN O'O'lex' are registered trademarks of Otometrics A/S. Echor-Screen', AUDX', Novigator' Pro, and ABaer' are registered trademarks of Natus Medical Incorporated.

Sanibel" is a registered trademark of Sanibel Supply'. All other Instrument brands or model names are copyrighted or trademarked by their respective owners and are used only for reference.

