

## MAICO Audiometer MA 25

### Light and independent for basic screening

#### Sleek and Small Design

The super-lightweight MA 25 weighs in at only 980 g. Small enough to hold in the palm of your hand! Tucked into its carrying case, the MA 25 is one of our smallest, lightest portables.

#### Accurate and Dependable Engineering

Because it is portable, the MA 25 has been engineered to be tough and reliable. Take the MA 25 everywhere with confidence. User-friendly knobs make adjusting frequency and hearing levels easy and comfortable!

#### Full Range of Test Options

The MA 25 offers a full range of testing opportunities. Frequencies range from 125 Hz to 8000 Hz and signal intensities run from  $-10 \text{ dB}_{\text{HL}}$  to  $100 \text{ dB}_{\text{HL}}$ . Sinus and warble tone (pulsed and continuous) make it complete.

#### Features at a glance

- Air Conduction with 11 frequencies from 125 to 8000 Hz
- Customize the frequencies to your own needs
- Can be operated with 3 standard AA batteries or eco-friendly with a medically approved external power supply (optional)
- Pure, pulse, warble tones
- Exceptionally lightweight, only 980 g
- MA 25e extends functionality of MA 25



## Technical Data

### MA 25 SPECIFICATIONS

<b>Test Signals</b>	Sinus and warble tone (pulsed and continuous)
<b>Test Frequencies Hz</b>	125, 250, 500, 750, 1000, 1500, 2000, 3000, 4000, 6000, 8000
<b>Level Steps</b>	5 or 1 dB
<b>Hearing Level Range with Headset DD45</b>	-10 dB <sub>HL</sub> - 100 dB <sub>HL</sub> 70 dB <sub>HL</sub> at 0.125 kHz, 90 dB <sub>HL</sub> at 0.25 and 8 kHz
<b>Audiometer Class</b>	Class 4 according to EN 60 645-1

### GENERAL

<b>Dimensions</b>	22.5 cm x 18 cm x 5.5 cm
<b>Weight</b>	980 g
<b>Power Supply</b>	3 AA batteries (100 hours life time typically) or optional power supply

### STANDARD COMPONENTS

- Audiogram Pad
- Operators Manual
- DD45 Headset Assembly

### OPTIONAL ACCESSORIES

- Hughson-Westlake test
- Patient response switch
- Power supply
- Carrying case
- Headset Holmco 8103, DD65 and TDH 39

### STANDARDS

EN 60 601-1  
EN 60 645-1, class 4,  
according to medical device directive  
93/42/EEC



### MA 25E EXTENDS FUNCTIONALITY OF MA 25

The MA 25e extends the MA 25 functionalities with an automatic test according to Hughson-Westlake. Parameters can easily be adapted to your personal preferences. The Talk Forward function makes the MA 25e easy to work with, particularly in sound booth installations.

**Standard  
Components**



MA 25 device




Headset DD45




Audiogram pad


**Optional  
Accessories**




Headset Holmco  
8103



Headset DD65



Headset TDH 39



Patient response  
switch

### MAICO Diagnostics GmbH

Sickingenstr. 70-71 · 10553 Berlin · Germany  
Tel.: +49 30/70 71 46-50 · Fax: +49 30/70 71 46-99  
sales@maico.biz · www.maico.biz