













WAX FILTER GUIDE

Product families	Ear piece type	Filter system	
		 ProWax	 ProWax miniFit
Intent™ Sirius™ Platform Real™ Polaris R™ Platform	miniFit Detect Receivers (60, 85, 100) miniFit Receivers (60, 80, 85, 100) 		✓
More™ Own™ Zircon Polaris™ Platform	LiteTip, Micro Mould, Power Flex Mould and MicroShell Detect (100, 105) 	✓	
Opn S™ Ruby Velox S™ Platform	MicroShell (85, Polaris R, Polaris), MicroShell Detect (60, 85) 		✓
Jet Opn™ Siya Velox Platform Alta Nera Ria Inium Platform	ITE, ITC, CIC, IIC 	Inium platform	Polaris & Velox platform Inium platform: non-wireless IIC & CIC V2 75 models
Agil Rise platform	Intiga ⁱ 	✓	
Acto Rise platform	S and M Receivers 	✓	
Ino Rise platform	P Receiver 		WaxStop
Intiga Rise platform	Micro Mould 		WaxStop
	Power Mould 	✓	
	ITE, ITC, CIC 	✓	

ProWax and **ProWax miniFit** are new, nano coated wax filters for ear pieces, providing better protection for longer performance.

ProWax can be used with all Oticon products and ear pieces where **NoWax** filters were previously used.

ProWax and ProWax miniFit are trademarks of Oticon.

Performance of the hearing devices is dependent on individual circumstances and may not be appropriate for all conditions. Always read the label and follow the instructions.