

People First

People First is our promise to empower people to communicate freely, interact naturally and participate actively

**Quality hearing care
by Oticon**

Quality Hearing Care

Would you like to boost your hearing in more everyday situation with a quality solution that's built to last?

Hearing is essential for communication, interaction and participation - not only at home but also in situations where noise makes listening more difficult. When hearing becomes difficult, not only do you need more sound - you need sound that is comfortable to listen to and lets you focus on important conversations.

113213UK / 06.11 Printed on 100% recycled paper.



Hearing instruments have come a long way

Modern digital hearing instruments handle everyday situations more efficiently and with fewer side effects than previous analogue aids. Annoying whistling and frequent volume adjustments are issues of the past. Now, you can enjoy good sound quality in contemporary designs that are robust and easy to handle.




Oticon Get

Oticon Get hearing instruments incorporate the latest technology by Oticon. The RISE technology is your guarantee for clear and comfortable sound quality that lets you focus on what you want to hear - conversations and important everyday events.

Depending on your needs, your hearing care professional will help choose the most discreet and comfortable solution for you. Oticon Get comes in a variety of styles that are elegantly designed and coloured to match your skin and hair for added discreetness.



	Behind-the-Ear	Completely-in-Canal	In-the-Canal	In-the-Canal Power	In-the-Ear
Prioritisation of speech from persons in front of you	●	-	●	Option	●
Reduction of noise for added listening comfort	●	●	●	●	●
Multiple listening programs	●	-	Option	Option	Option
Volume control	●	-	-	Option	Option
 Teleloop compatible	●	-	Option	Option	Option
Audio input connector for attachment of FM receiver	●	-	-	-	-
Protective surface coating	●	-	-	-	-

Behind-the-Ear instrument colours



by Oticon

Get is not made by just anyone. Oticon is a world renowned manufacturer of high-performance hearing instruments for people who are seeking the very best.