





# Light and Easy

Looking for a light solution to boost everyday audibility for your clients? One that's fast and easy to fit? Oticon Get is it!

## **Hearing Care**

Oticon Get offers some of the most fundamental benefits of hearing care: audibility with open ear acoustics and directionality - plus comfort in noisy environments.

# by Oticon

Where business relationships are supportive and long lasting.

> Get is Oticon's most affordable entry to proven and reliable BTEs and custom solutions.

### **Main Features**

- Open Ear Acoustics
- Corda2 thin tube and hook versions
- Nano-coated BTEs
- Surround, Split and Full Directionality (manual)
- Noise Management (modulation)
- Dynamic Feedback Cancellation DFC2
- 4 user programs
- Telecoil, DAI & FM options

### Fitting flexibility

- Fitting bandwidth 6.5 kHz
- 4 fitting bands
- NAL-NL1, NAL-NL2, DSL v5.0a m[i/o]
- Adaptation Manager
- In-situ Audiometry

#### CUSTOM





ITC Power (Omni or Dir)

BTE 13




Corda<sup>2</sup> thin tube












Earmould Corda<sup>2</sup>



ITC Power





BTE POWER







Chroma Beige

Diamond Black

Steel Grey

Chestnut Brown Terracotta

## People First

People First is our promise to empower people to communicate freely, interact naturally and participate actively

