

Oticon charger compatibility overview



Oticon SmartCharger miniRITE Oticon Charger miniRITE Oticon SmartCharger miniRITE R Oticon SmartCharger miniBTE R Charger 1.0 miniRITE R Charger 1.0 miniBTE R

Oticon Intent™	•	•	-	-	-	-
Oticon Real™	-	-	•	•	•	•
Oticon More™	-	-	•	•	•	•
Oticon Zircon	-	-	•	•	•	•
Oticon Opn S™	-	-	-	-	•	-
Oticon Ruby	-	-	-	-	•	-
Oticon Play PX	-	-	•	•	•	•
Oticon Opn Play™	-	-	-	-	•	-
Oticon CROS PX	-	-	•	-	•	-

Performance of the hearing devices is dependent on individual circumstances and may not be appropriate for all conditions. Always read the label and follow the instructions.