

VIDEO  
**OTOSCOPE**

The GSI Video Otoscope is intended to allow inspection of the external ear canal and tympanic membrane under magnification. These examinations are useful in the evaluation, identification, documentation, and diagnosis of ear disorders. It is also used to assist with hearing aid and probe microphone placement for Real Ear Measurements.

**KEY FEATURES**

- Lightweight
- USB-powered
- Convenient capture button
- Generate clear images of the ear canal and tympanic membrane
- Images immediately available for patient counselling
- Versatile functionality by operating directly via GSI Suite™
- Images can be transferred to Noah with GSI Suite module

